

Phonak

Fast Facts

Feedback threshold overtuning





What is it?

Feedback threshold overtuning provides access to additional high frequency gain without having to run a feedback test. Traditionally this gain has been limited by the predicted feedback threshold.

It is highly recommended that a feedback measurement is completed as standard practice however, overtuning can be beneficial in the following situations:

- When a feedback measurement is not possible; or
- If access to additional high frequency gain is needed after running the measurement

How to use it?

<p>Step 1. Run a feedback measurement</p> 	<p>Step 2. Optimize acoustics and repeat measurement if needed based on feedback measurement</p> <ul style="list-style-type: none"> • Overtuning does not replace managing feedback with a more occluding earpiece 
<p>Step 3. Adjust overtuning in Feedback & real ear test if needed</p> 	<p>Step 4. Check sound quality</p> 

Important notes for Step 3:

- Even though an additional 12 dB of gain is possible, it is recommended that overtuning is kept to 5 dB or less.
 - This is to manage feedback risk, sound quality and battery consumption
 - A warning will appear in Phonak Target when 5 dB is exceeded
- The number of overtuning steps available is dependent on the hearing loss and the acoustics parameters
- The number of overtuning steps does not indicate the additional gain in dB already applied