

Phonak Audéo™ P

Product information



Phonak Audéo P is the comprehensive RIC portfolio for mild to profound hearing losses. Phonak Audéo P includes four designs, four performance levels and five external receivers, including the new ActiveVent™ receiver.

Powered by the innovative technology of the Phonak Paradise platform, Phonak Audéo P hearing aids are the solution for those clients that prefer flexible and discreet RIC housings.

For more information visit www.phonakpro.com.au

Product description



Audéo P-312
Audéo P-312 Trial
P90/P70/P50/P30



Audéo P-13T
Audéo P-13T Trial
P90/P70/P50/P30



Audéo P-R
Audéo P-R Trial
P90/P70/P50/P30



Audéo P-RT
Audéo P-RT Trial
P90/P70/P50/P30

Battery size		312	13		
Multi-function button		•	•	•	•
Telecoil			•		•
Bluetooth®		•	•	•	•
Nano coating		•	•	•	•
IP rating		IP68 ¹	IP68 ¹	IP68 ¹	IP68 ¹
Dimensions	L x W x D	26.8 x 11.2 x 7.1 mm	29.0 x 11.9 x 8.1 mm	28.0 x 11.9 x 8.0 mm	30.5 x 12.5 x 8.0 mm
	Weight	1.1 g	1.4 g	2.1 g	2.4 g
Charging times ²	0% to 100%			3 h	3 h
	0% to 80%			1.5 h	1.5 h
	30 min. fitting			10 min.	10 min.

S Receiver

Max. Power Output (dB SPL)	2 cc coupler	111	111	111	111
Max. gain (dB)	2 cc coupler	47	47	47	47
Frequency range (Hz)		<100 - >8000	<100 - >8000	<100 - >8000	<100 - >8000
Working current (mA)		1.4	1.4		

M Receiver

Maximum Power Output (dB SPL)	2 cc coupler	114	114	114	114
Max. gain (dB)	2 cc coupler	51	51	51	51
Frequency range (Hz)		<100 - >8000	<100 - >8000	<100 - >8000	<100 - >8000
Working current (mA)		1.4	1.4		

MAV Receiver

Maximum Power Output (dB SPL)	2 cc coupler			119	119
Max. gain (dB)	2 cc coupler			50	50
Frequency range (Hz)				<100 - >8000	<100 - >8000
Working current (mA)					

P Receiver

Maximum Power Output (dB SPL)	2 cc coupler	122	122	122	122
Max. gain (dB)	2 cc coupler	59	59	59	59
Frequency range (Hz)		<100 - 6300	<100 - 6300	<100 - 6300	<100 - 6300
Working current (mA)		1.5	1.5		

UP Receiver

Maximum Power Output (dB SPL)	2 cc coupler	132	132	132	132
Max. gain (dB)	2 cc coupler	71	71	71	71
Frequency range (Hz)		<100 - 6100	<100 - 6100	<100 - 6100	<100 - 6100
Working current (mA)		1.4	1.4		

¹ IP68 indicates that the hearing aid is water and dust resistant. It survived continuous immersion in 1 meter of fresh water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard.

² Actual charging time depends on the remaining battery life, but it will not exceed 3 hours.

Housing colours



P1
Sand Beige



H0
Beige



P3
Sandalwood



P4
Chestnut



P5
Champagne



P6
Silver Gray



P7
Graphite Gray



P8
Velvet Black

Performance levels

	Premium	Advanced	Standard	Essential
AutoSense OS™ 4.0	Premium	Advanced	Standard	Essential
Comfort in Echo	•			
Speech in Car	•			
Speech in Loud Noise	•			
Music	•	•		
Comfort in Noise	•	•	•	
Speech in Noise	•	•	•	•
Calm Situation	•	•	•	•
Media speech + mic	•	•	•	•
Media music + mic	•	•	•	•
Additional programs				
Max. additional programs	4	4	4	4
Comfort in Echo	•			
Speech in Loud Noise	•	•		
Speech in 360°	•	•		
Speech in Noise	•	•	•	•
Calm Situation	•	•	•	•
Music	•	•	•	•
Acoustic phone	•	•	•	•
Custom program	•	•	•	•
Streaming programs				
RogerDirect + mic	•	•	•	•
PartnerMic + mic	•	•	•	•
Phone call + mic	•	•	•	•

Features	Premium	Advanced	Standard	Essential
UltraZoom	•	•	•	•
Dynamic Noise Cancellation	•			
Speech Enhancer	•			
EchoBlock	•			
Tap Control ^{a)} - Siri/Google Assistant	•	•		
- Pause/resume streaming	•	•		
- Accept/end phone call	•	•		
StereoZoom	•	•		
WindBlock	•	•		
SoundRelax	•	•		
Motion Sensor Hearing ^{a)}	•	•	•	
DuoPhone	•	•	•	
Real Ear Sound	•	•	•	
Roger and Directional setting	•	•	•	
SoundRecover2	•	•	•	•
WhistleBlock	•	•	•	•
NoiseBlock	•	•	•	•
RogerDirect™	•	•	•	•
User Preference Tuning	•	•	•	•
QuickSync	•	•	•	•
AOV	•	•	•	•
Tinnitus Balance	•	•	•	•
auto Acclimatization	•	•	•	•
Environmental balance	•	•	•	•
Fine tuning channels	20	20	16	12

^{a)} Only available for Phonak Audéo P-R and RT

Fitting

	Audéo P-312 Audéo P-312 Trial	Audéo P-13T Audéo P-13T Trial	Audéo P-R Audéo P-R Trial	Audéo P-RT Audéo P-RT Trial
Phonak Target	7.0 or higher	7.0 or higher	7.1.9 or higher	7.1.9 or higher
Noahlink Wireless	•	•	•	•
NOAHLink		•		
HI-PRO, HI-PRO2		•		

Acoustic coupling

Length	S Receiver	M Receiver	MAV Receiver	P Receiver	UP Receiver
Available in four different lengths (0, 1, 2, 3), left and right Dome	•	•	•	•	•
Dome					
Cap dome	•	•		•	
Open dome Available in three sizes (S, M, L)	•	•		•	
Vented dome Available in three sizes (S, M, L)	•	•		•	
Power dome Available in three sizes (S, M, L)	•	•		•	
SlimTip*					
Acryl	•	•		•	
Silicone	•	•		•	
Titanium AV			•		
Universal AV Available in four shapes (A, B, C, D)			•		
cShell					
Acryl	•	•		•	•
Titanium	•	•		•	•
Fitting ranges					

Mild to moderately-severe hearing loss.

Mild to severe hearing loss.

Mild to severe hearing loss.

Mild to severe hearing loss.

Moderate to profound hearing loss.

*Available with extra canal length option

ActiveVent™ Receiver

The brand new ActiveVent Receiver (MAV) offers the benefit of both open and closed fit in one receiver. It is compatible with all rechargeable Audéo P devices. ActiveVent receivers can be instantly fit for in-clinic demonstration or short term home trial, with a Universal SlimTip AV that is available in 4 generic shapes. A SlimTip Titanium AV that is custom made for the client's ear is recommended for long term wear and to achieve optimal acoustic performance.

Phonak Audéo P Trial™

The Phonak Audéo P Trial hearing aids can be quickly set to the required performance level in order to suit your recommendation. All levels of technology in one Phonak Audéo P Trial hearing aid provides you with many options for providing the hearing solution at the first appointment or bridging the gap during repairs.

Phonak Audéo P charging options



Phonak Charger Case Combi

The Phonak Charger Case Combi is a charger, drying kit and protective hard case all in one, offering a robust and splash-proof housing.



Phonak Power Pack

For charging rechargeable hearing aids on the go, the Phonak Power Pack can be easily attached to the bottom of the Phonak Charger Case Combi.



Phonak Mini Charger Case

The Phonak Mini Charger Case, with standard USB-C charging outlet, is a compact charger for two hearing aids.

Phonak Audéo P Bluetooth® versions

Direct calls/audio streaming	Bluetooth® 4.2 wireless technology and most older Bluetooth phones
------------------------------	--

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Sonova AG is under license.

Wireless communication portfolio

	Audéo P-312 Audéo P-312 Trial P90/P70/P50/P30	Audéo P-13T Audéo P-13T Trial P90/P70/P50/P30	Audéo P-R Audéo P-R Trial P90/P70/P50/P30	Audéo P-RT Audéo P-RT Trial P90/P70/P50/P30
Wireless accessory				
TV Connector	•	•	•	•
Phonak PartnerMic	•	•	•	•
Phonak RemoteControl	•	•	•	•
Roger				
RogerDirect™	•	•	•	•
Apps				
myPhonak	•	•	•	•

For more information please visit www.phonakpro.com/audeoparadise
For information on Roger microphones please visit www.phonakpro.com/roger
For more information on eSolutions please visit www.phonakpro.com/smartapps

life is on

At Phonak, we believe that well-hearing equates to well-being and is essential to living life to the fullest. For more than 70 years, we have remained passionate about creating a world where 'life is on' for everyone. Our innovative hearing solutions are designed for people of all ages and all degrees of hearing loss, to connect socially, thrive mentally and emotionally.

www.phonak.com.au



027-0602-02/V3.00/2021-04/na © 2021 Sonova AG All rights reserved

