

Your Roger On Sphere™ Welcome Guide.

Congratulations on starting your Roger On Sphere trial!

With Roger on your side, you're equipped with a powerful and versatile hearing solution that will upgrade your listening experience and empower you to be confident in more places.

This guide will walk you through everything you need to know to get the most from your new Roger microphone. Get ready to experience and expand your world with Roger On Sphere!

Getting started with your Roger On

1 Set up RogerDirect™

This is a one-time step to connect Roger with your hearing aids. Your audiologist may have already done this during your last clinic visit.

If not, refer to the step-by-step instructions in the "Roger receiver setup" guide included in the box or watch the Roger "How-To" video on the Phonak YouTube channel

2 Get familiar with Roger On

As you unbox your new microphone, take a moment to go through the Quick setup guide provided in the box. It covers what's inside, explains the icons that appear on the screen, and gives you examples of everyday usage.

3 Charge before use

It takes 2 hours to fully charge. A full battery can provide up to 8 hours of continuous use.

4 Download the myRogerMic app

The app lets you adjust Roger's settings easily. It's optional but a great way to enhance your experience.



Google Play



App Store

And that's it - you're all set and ready to bring Roger along wherever you go.



roger

Put Roger to the test

Try it in different settings and note down whether it helped (Yes/ No). When your trial period ends, we recommend sharing your findings with your audiologist to get the most from your experience.

At social gatherings	Recommended use	Yes	No
At the coffee shop/ a noisy restaurant, dining with friends and family	Place Roger On Sphere in the centre of the table. You can select the microphone direction	☐	☐
At a social event, with loud music and people talking	Hold it in your hand and point it towards the person speaking	☐	☐
At the shopping mall/ grocery store	Hold it in your hand and point it towards the person speaking	☐	☐
Sports and other activities			
At the gym or training outside	Give Roger On Sphere to your coach/ teammate and ask them to clip it on their shirt. Alternatively, they can wear it around their neck by using the lanyard	☐	☐
Hiking or biking	Give it to a fellow hiker/ biker to clip on their shirt while you're on the go	☐	☐
At home			
Conversation with your partner or a friend when you're not face-to-face (e.g. in a different room)	Give Roger On Sphere to your partner/ friend to clip it on their shirt	☐	☐
Watching TV	Use the docking station to connect Roger to your TV	☐	☐
When travelling			
Have a conversation with your travel companion, while in the car or in public transportation	Give Roger On Sphere to your travel companion to clip it on their shirt	☐	☐
Business events			
Lectures and presentations	Give Roger On to the presenter and sit comfortably among the audience	☐	☐
Panel discussion	Place it on the stage, ideally on a table, where the panel discussion takes place. Open the mobile app and choose the Wide pointing mode. You can now sit comfortably among the audience and listen to the discussion	☐	☐
At work			
Team workshops	Place Roger On Sphere in the middle of the table when having a team discussion. Give it to the instructor to clip it on their shirt. You can also connect multiple Roger microphones together, through the mobile app, with the feature "MultiTalker Network"	☐	☐
Standing meetings and hallway conversations	Hold it in your hand and point it towards the person speaking	☐	☐
Small meetings	Place it in the middle of the meeting room table	☐	☐
Online calls	You can also use Roger as a headset. First, use the cable to connect it to your laptop. Then clip it on your shirt. It will pick up your voice and send it to the computer	☐	☐

