

Phonak CROS II Order Form

Account No: 200 _____









Hearing Clinic:	Client Name:
Delivery Address:	Date Ordered:
Contact Name:	Date Required:
Contact Number:	VIP: 24hr in house \$50 per order (+GST)
Email:	

Step 1: Choose your Phonak CROS II transmitter for the unaidable ear

Behind-The-Ear Options	Order No.	OHS Code	Colour*	Qty
CROS II-312	050-0228-xx	N0091AD	_____	_____
CROS II-13	050-0229-xx	N0092AD	_____	_____

*See overleaf for colour options

Step 2: Choose your Acoustic Coupling Option (For more information, refer to 'Earpeice Options' on following page)

Phonak CROS Demo Hook* <input checked="" type="checkbox"/> L <input checked="" type="checkbox"/> R  <small>*CROS Demo Hook included as standard</small>	CROS Slim Tube with Domes  L R  Indicate Tube Length 0-3	CROS Slim Tube with Slim Tip^  L R  Indicate Tube Length 0-3 ^ Impression Required	CROS Tip^  L R  Indicate Tube Length 0-3 ^ Impression Required	CROS B Mould^1  L R ^ Impression Required ^1 Compatible with CROS II-13 Only. Comes with CROS II Hook
---	---	--	---	---

Step 3: Choose from the following options

CROS Slim Tubes	Open Domes	Closed Domes	Power Domes	Colour & Options
0-L 054-0690	0-R 054-0686	S 054-1987	S 054-1990	S 054-1993
1-L 054-0691	1-R 054-0687	M 054-1988	M 054-1991	M 054-1994
2-L 054-0692	2-R 054-0688	L 054-1989	L 054-1992	L 054-1995
3-L 054-0693	3-R 054-0689			

Colour & Options: Transparent (21) Removal Line
Pink (26) Extra canal length needed?
S = Standard

Step 4: Choose your Wireless Venture hearing instrument for the better ear

Audéo V Please complete the Phonak Audéo V RIC Order Form to specify correct inclusions.	Bolero V Please complete the Phonak Bolero V BTE Order Form to specify correct inclusions.
--	--



Step 5: Choose your Wireless Communication Products

Refer to Phonak's Wireless Communication Order Form

Step 6: Check List







- RIC Order Form accompanies this order
- BTE Order Form accompanies this order
- Wireless Communication Order Form accompanies this order

Product description

			
	CROS II-312	CROS II-13	
Battery size	312	13	¹ IP57 indicates that the hearing aid is water resistant and dust protected. It survived immersion in 1m of water for 30 minutes and 8 hours in a dust chamber as per the IEC60529 standard.
Push button	•	•	
Volume control	•	•	² IP67 indicates that the hearing aid is water resistant and dust tight. It survived immersion in 1m of water for 30 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.
Nano coating	•	•	
IP rating	IP57 ¹	IP67 ²	
Dimensions (L x W x D)	2.4 x 0.6 x 0.9 cm (0.94 x 0.24 x 0.35")	3.2 x 0.8 x 1.5 cm (1.26 x 0.31 x 0.59")	
Weight	1 g (0.035 oz)	2.2 g (0.078 oz)	

Housing colours

Hair & Skin

			
P1	P2	P3	P4
Sand Beige	Amber Beige	Sandalwood	Chestnut
			
P5	P6	P7	P8
Champagne	Silver Gray	Graphite Gray	Velvet Black

Fashion

	
P9	Q1
Ruby	Petrol

Traditional


O1
Beige

Earpiece options

CROS DemoHook	<ul style="list-style-type: none">• Hook to hold CROS on the ear during fitting and demonstration• For demo purposes only• Included with any Phonak CROS and CROS II device replacing the CROS retentions
CROS SlimTube	<ul style="list-style-type: none">• Standard SlimTube for CROS devices• Available in left/right and four different lengths (0, 1, 2, 3)• Domes need to be attached. Available as open, closed and power dome in three sizes (S, M, L)
CROS SlimTip	<ul style="list-style-type: none">• Custom made earpiece for CROS devices• Available in left/right and four different lengths (0, 1, 2, 3)• Available in hard and soft material
CROS II Hook	<ul style="list-style-type: none">• Earhook to attach a custom earmold with standard tubing to CROS II-13• Only available for the CROS II-13 model• Included with any CROS II-13 device or can be ordered separately (order number: 048-0261)