



Phonak Lyric™

The contact lens for your ear.™

Lyric is the world's only 100% invisible, extended-wear hearing device you can wear 24/7 for months at a time.¹



Actual Size

PHONAK
life is on



Phonak Lyric™

"I completely forget about my hearing loss because you never take it out, you never put it back, you feel like a regular person. You're all whole again."

– Pete, Lyric Wearer

93% of Lyric users would recommend Lyric to a friend or loved one.²

100% invisible

Revolutionary design allows Lyric to sit invisibly in the ear canal, out of sight and out of mind.

Clear, natural sound

Captures sound naturally even in noisy background situations so you hear everything that counts.


No daily hassles

Designed to be worn 24/7 for months at a time.¹

No batteries to change

No daily insertion, removal, or battery changes required. Ever.






How Lyric works

Lyric is placed comfortably close to the ear drum. The placement of Lyric enables the ear to capture sound and funnel it toward the ear drum organically, just as it would without a hearing device. As a result, Lyric is able to offer pure, natural sound quality and improved hearing in noisy background situations.



www.phonak.com/lyric



"I was able to hear sounds I hadn't heard in years... It was just like turning on a switch. It was amazing."

– Libby,
Lyric Wearer

86% of Lyric users agree
Lyric's sound quality is very natural.³

Benefits of Lyric

Unlike other hearing solutions, which are inserted and removed on a daily basis, extended-wear Lyric hearing devices are worn around the clock for months at a time.¹ Users can wear Lyric during most daily activities, like sleeping, showering, and exercising⁴ And, Lyric never requires battery changes, cleaning or repairs.

www.phonak.com/lyric

Contact your local authorised Lyric provider:

¹Individual replacement needs may vary. Duration of device battery life varies by patient and is subject to individual ear conditions. ²Based on results from two surveys totaling 109 patients who had worn Lyric for at least 30 days. ³Based on results from a survey totaling 113 patients who had worn Lyric for at least 30 days. ⁴Lyric is water resistant, not waterproof, and should not be completely submerged under water.

028-1175-02. ©Phonak AG. All rights reserved.

