



Phonak Product Catalogue

Effective February 2022

A Sonova brand

PHONAK
life is on

Table of Content

Hearing Instruments:

Paradise Platform 3

Phonak Audéo™ P	4
Phonak Audéo™ P Rechargeable	5
Phonak CROS P	6
ActiveVent™ Receiver	7
Phonak Naída™ P	8
Phonak Chargers: Compatibility Overview	9

Hearing Instruments:

Marvel Platform 10

Phonak Audéo™ M	11
Phonak Audéo™ M Rechargeable	12
Phonak Bolero™ M	13
Phonak Bolero™ M Rechargeable	14
Phonak Naída™ M	15
Phonak Virto™ M	16
Phonak Virto™ M-Titanium	17
Phonak Sky™ M	18

Hearing Instruments:

Belong Platform 19

Phonak Audéo™ B	20
Phonak Audéo™ B Rechargeable	21
Phonak Bolero™ B	22
Phonak Bolero™ B Rechargeable	23
Phonak CROS B	24
Phonak CROS B Rechargeable	25
Phonak Naída™ B	26
Phonak Naída™ B Rechargeable	27
Phonak Naída™ Link	28
Phonak Sky™ B	29
Phonak Sky™ B Rechargeable	30
Phonak Virto™ B	31

Hearing Instruments:

Free-to-Client 32

Phonak Vitus	33
--------------------	----

Hearing Protection 34

Phonak Serenity Choice™	35
Phonak Serenity Choice™ Plus	36

Programming Tools 37

Programming Tools	38
-------------------------	----

Wireless Accessories 39

Phonak Wireless Communication Portfolio	40
Roger™ for Marvel and Paradise	42
Roger™ Portfolio	43
Roger™ Portfolio for Education	45
Roger™ Dynamic Soundfield	46
Phonak Colour Chart	48

Sound Delivery Systems 49

SDS 4.0 Receivers & SlimTubes	50
SDS 4.0 Domes & Custom Coupling	51
Receivers, Domes & Coupling	52
SlimTubes	54
Tone Hooks	55
Cleaning and Wax Management	56
Microphone Protection	57
Batteries	58
Consumables	59
Tools	60

Policies and purchase conditions 61

Policies and purchase conditions	62
--	----

Hearing Instruments: Paradise Platform

Phonak Audéo™ P

Premium (90), Advanced (70), Standard (50), Essential (30)



Audéo P-312



Audéo P-13T

Model	Order Reference Refer to colour chart	HSP Code	Category	Gain / MPO (2cc Coupler)				Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS P Compatibility
				S	M	P	UP						
Audéo P90-312*	050-0769-xx	B679PH0	C2	47/111	51/114	59/122	71/132	312	IP68 ¹	Y	N	Y	N
Audéo P70-312*	050-0788-xx	B680PH0	C2	47/111	51/114	59/122	71/132	312	IP68 ¹	Y	N	Y	N
Audéo P50-312*	050-0787-xx	B681PH0	C2	47/111	51/114	59/122	71/132	312	IP68 ¹	Y	N	Y	N
Audéo P30-312*	050-0786-xx	B682PH0	C2	47/111	51/114	59/122	71/132	312	IP68 ¹	Y	N	Y	N
Audéo P90-13T	050-0770-xx	B683PH0	C2	47/111	51/114	59/122	71/132	13	IP68 ¹	Y	Y	Y	Y
Audéo P70-13T	050-0791-xx	B684PH0	C2	47/111	51/114	59/122	71/132	13	IP68 ¹	Y	Y	Y	Y
Audéo P50-13T	050-0790-xx	B685PH0	C2	47/111	51/114	59/122	71/132	13	IP68 ¹	Y	Y	Y	Y
Audéo P30-13T	050-0789-xx	B686PH0	C2	47/111	51/114	59/122	71/132	13	IP68 ¹	Y	Y	Y	Y
Audéo P-312 Trial	050-0844-P6	N/A	N/A	47/111	51/114	59/122	71/132	312	IP68 ¹	Y	N	Y	N
Audéo P-13T Trial	050-0845-P6	N/A	N/A	47/111	51/114	59/122	71/132	13	IP68 ¹	Y	Y	Y	Y

* Noalink wireless essential for programming.

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Colour chart



P1

Sand Beige



P3

Sandalwood



P4

Chestnut



P5

Champagne



P6

Silver Grey



P7

Graphite Grey



P8

Velvet Black



H0

Beige

Phonak Audéo™ P Rechargeable

Premium (90), Advanced (70), Standard (50), Essential (30)



Audéo P-R



Audéo P-RT

Model	Order Reference Refer to colour chart	HSP Code	Category	Gain / MPO (2cc Coupler)				Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS P Compatibility
				S	M	P	UP						
Audéo P90-R*	050-0771-xx	B687PHO	C2	47/111	51/114	59/122	71/132	Rechargeable	IP68 ¹	Y	N	Y	Y
Audéo P70-R*	050-0794-xx	B688PHO	C2	47/111	51/114	59/122	71/132	Rechargeable	IP68 ¹	Y	N	Y	Y
Audéo P50-R*	050-0793-xx	B689PHO	C2	47/111	51/114	59/122	71/132	Rechargeable	IP68 ¹	Y	N	Y	Y
Audéo P30-R*	050-0792-xx	B690PHO	C2	47/111	51/114	59/122	71/132	Rechargeable	IP68 ¹	Y	N	Y	Y
Audéo P90-RT*	050-0772-xx	B691PHO	C2	47/111	51/114	59/122	71/132	Rechargeable	IP68 ¹	Y	Y	Y	Y
Audéo P70-RT*	050-0797-xx	B692PHO	C2	47/111	51/114	59/122	71/132	Rechargeable	IP68 ¹	Y	Y	Y	Y
Audéo P50-RT*	050-0796-xx	B693PHO	C2	47/111	51/114	59/122	71/132	Rechargeable	IP68 ¹	Y	Y	Y	Y
Audéo P30-RT*	050-0795-xx	B694PHO	C2	47/111	51/114	59/122	71/132	Rechargeable	IP68 ¹	Y	Y	Y	Y
Audéo P-R Trial	050-0846-P6	N/A	N/A	47/111	51/114	59/122	71/132	Rechargeable	IP68 ¹	Y	N	Y	Y
Audéo P-RT Trial	050-0847-P6	N/A	N/A	47/111	51/114	59/122	71/132	Rechargeable	IP68 ¹	Y	Y	Y	Y

* Noalink wireless essential for programming.

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Charging options

Phonak Charger BTE RIC 075-3016-10

This charger is suitable for all Paradise rechargeable BTE and RIC hearing aids. It features a holding space for large earmoulds, a USB-C charging inlet and an indicator light.



Phonak Charge and Care 075-0041-10

The Phonak Charge and Care is a charger that dries hearing aids with active ventilated air flow and treats hearing aids with UV light.



Phonak Charger Case Combi 075-0040

The Phonak Charger Case Combi is a charger, drying kit and protective hard case all in one.



Phonak Power Pack 071-0003

For charging rechargeable hearing aids on the go, the Phonak Power Pack can be easily attached to the bottom of the Phonak Charger Case Combi.



Colour chart

P1	P3	P4	P5	P6	P7	P8	H0
Sand Beige	Sandalwood	Chestnut	Champagne	Silver Grey	Graphite Grey	Velvet Black	Beige

Phonak CROS P

Compatible with Audéo™ Paradise



Audéo P-R

Audéo P-RT

Model	Order Reference Refer to colour chart	HSP Code	Category	Battery Size	IP Rating	Compatibility
CROS P-13	050-0861-xx	N0275CR	NS	13	IP68 ¹	Audéo P-13T
CROS P-R	050-0860-xx	N0276CR	NS	Rechargeable	IP68 ¹	Audéo P-RT and Audéo P-R

* Noahlink wireless essential for programming.

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

CROS P performance

CROS P offers the Speech in Loud Noise program with StereoZoom. It is available in combination with hearing aids of Premium and Advanced performance levels. For the Standard and Essential performance levels, CROS P provides Real Ear Sound and omnidirectional microphone modes. On the better ear, CROS P clients fully benefit from all Paradise features.

Charging options

Phonak Charger Case Combi 075-0040

The Phonak Charger Case Combi is a charger, drying kit and protective hard case all in one.



Phonak Power Pack 071-0003

For charging rechargeable hearing aids on the go, the Phonak Power Pack can be easily attached to the bottom of the Phonak Charger Case Combi.



Colour chart

P1 Sand Beige	P3 Sandalwood	P4 Chestnut	P5 Champagne
P6 Silver Grey	P7 Graphite Grey	P8 Velvet Black	H0 Beige

	Order Reference
CROS Tube 4.0 0L	054-0886
CROS Tube 4.0 0R	054-0887
CROS Tube 4.0 1L	054-0888
CROS Tube 4.0 1R	054-0889
CROS Tube 4.0 2L	054-0890
CROS Tube 4.0 2R	054-0891
CROS Tube 4.0 3L	054-0892
CROS Tube 4.0 3R	054-0893
SlimTip 4.0 Acrylic	062-0006-01
SlimTip 4.0 Silicone	062-0007-01
CROS Tip 4.0 Acrylic	062-0016-**

Earpiece options

Standard dome
(CROS Slim Tube 4.0)

- Standard dome available in three sizes (S, M, L)
- CROS Slim Tube 4.0 is available in left/right and four different lengths (0, 1, 2, 3)

SlimTip (CROS Slim Tube 4.0)

- Custom made earpiece
- CROS Slim Tube 4.0 is available in left/right and four different lengths (0, 1, 2, 3)

CROS Tip 4.0

- Fully custom made earpiece for CROS, available in four different cable lengths (0, 1, 2, 3)

ActiveVent™ Receiver



ActiveVent Receivers	Order Reference
MAV Receiver 0L	054-0962
MAV Receiver 1L	054-0963
MAV Receiver 2L	054-0964
MAV Receiver 3L	054-0965
MAV Receiver 0R	054-0966
MAV Receiver 1R	054-0967
MAV Receiver 2R	054-0968
MAV Receiver 3R	054-0969
Fit'n'Go Kit	
ActiveVent Fit'n'Go Kit	054-0991

Compatible with Audéo Paradise rechargeable hearing aids (requires latest firmware.)
ActiveVent Receiver comes with a 6 month warranty.

Titanium Slim Tip	Order Reference
SlimTip Titanium AV	062-0049-01
Universal SlimTip	
Universal SlimTip AV A-L	054-0998
Universal SlimTip AV A-R	054-1002
Universal SlimTip AV B-L	054-0997
Universal SlimTip AV B-R	054-1001
Universal SlimTip AV C-L	054-0996
Universal SlimTip AV C-R	054-1000
Universal SlimTip AV D-L	054-0995
Universal SlimTip AV D-R	054-0999
Wax filter	
Filter HF3 red (Dispenser)	098-0292
Filter HF3 blue (Dispenser)	098-0293

Phonak Naída™ P

Premium (90), Advanced (70), Standard (50), Essential (30)



Naída P-PR



Naída P-UP

Model	Order Reference Refer to colour chart	HSP Code	Category	Gain / MPO (2cc Coupler)		Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS P Compatibility
				SlimTube	Earhook						
Naída P90-PR*	050-0810-xx	B711PH0	C2	69/129	68/130	Rechargeable	IP68 ¹	Y	N	Y	N
Naída P70-PR*	050-0809-xx	B710PH0	C2	69/129	68/130	Rechargeable	IP68 ¹	Y	N	Y	N
Naída P50-PR*	050-0808-xx	B709PH0	C2	69/129	68/130	Rechargeable	IP68 ¹	Y	N	Y	N
Naída P30-PR*	050-0807-xx	B708PH0	C2	69/129	68/130	Rechargeable	IP68 ¹	Y	N	Y	N
Naída P90-UP*	050-0818-xx	B715PH0	C1	N/A	84/141	675	IP68 ¹	Y	Y	Y	N
Naída P70-UP*	050-0817-xx	B714PH0	C1	N/A	84/141	675	IP68 ¹	Y	Y	Y	N
Naída P50-UP*	050-0816-xx	B713PH0	C1	N/A	84/141	675	IP68 ¹	Y	Y	Y	N
Naída P30-UP*	050-0815-xx	B712PH0	C1	N/A	84/141	675	IP68 ¹	Y	Y	Y	N
Naída P-PR Trial	050-0837-P6	N/A	N/A	69/129	68/130	Rechargeable	IP68 ¹	Y	N	Y	N
Naída P-UP Trial	050-0839-P6	N/A	N/A	N/A	84/141	675	IP68 ¹	Y	Y	Y	N

* Noalink wireless essential for programming.

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Charging option

Phonak Charger BTE RIC
075-3016-10



This charger is suitable for all Paradise rechargeable BTE and RIC hearing aids. It features a holding space for large earmoulds, a USB-C charging inlet and an indicator light.

Colour chart

P1	P3	P4	P5	P6	P7	P8	H0
Sand Beige	Sandalwood	Chestnut	Champagne	Silver Grey	Graphite Grey	Velvet Black	Beige

Phonak Chargers: Compatibility Overview

	Phonak Charge and Care	Phonak Charger BTE RIC	Phonak Charger Case Combi	Phonak Power Pack
				
	075-0041-xx	075-3016-xx	075-0040-xx	071-0003
Phonak Audéo P-R	•	•		•
Phonak Audéo P-RT	•	•		•
Phonak CROS P-R		•		•
Phonak Naída P-PR		•		
Phonak Audéo M-R	•	•		•
Phonak Audéo M-RT	•	•		•
Phonak Sky M-PR		•		
Phonak Bolero M-PR		•		
Phonak Audéo B-R				•
Phonak Bolero B-PR				•
Phonak Sky B-PR				•
Phonak Naída B-R RIC				•
Phonak CROS B-R				•
Roger Focus II		•		

Hearing Instruments: Marvel Platform

Phonak Audéo™ M

Premium (90), Advanced (70), Standard (50), Essential (30)



Audéo M-312



Audéo M-13T



Audéo M-312T

Model	Order Reference Refer to colour chart	HSP Code	Category	Gain / MPO (2cc Coupler)				Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS B Compatibility
				S	M	P	UP						
Audéo M90-312*	050-0494-xx	B599PHO	C2	46/111	50/114	58/112	67/130	312	IP68 ¹	Y	N	Y	N
Audéo M70-312*	050-0493-xx	B601PHO	C2	46/111	50/114	58/112	67/130	312	IP68 ¹	Y	N	Y	N
Audéo M50-312*	050-0492-xx	B603PHO	C2	46/111	50/114	58/112	67/130	312	IP68 ¹	Y	N	Y	N
Audéo M30-312*	050-0491-xx	B605PHO	C2	46/111	50/114	58/112	67/130	312	IP68 ¹	Y	N	Y	N
Audéo M90-13T	050-0519-xx	B615PHO	C2	46/111	50/114	58/112	67/130	13	IP68 ¹	Y	Y	Y	N
Audéo M70-13T	050-0518-xx	B616PHO	C2	46/111	50/114	58/112	67/130	13	IP68 ¹	Y	Y	Y	N
Audéo M50-13T	050-0517-xx	B617PHO	C2	46/111	50/114	58/112	67/130	13	IP68 ¹	Y	Y	Y	N
Audéo M30-13T	050-0516-xx	B618PHO	C2	46/111	50/114	58/112	67/130	13	IP68 ¹	Y	Y	Y	N
Audéo M90-312T	050-0498-xx	B611PHO	C2	46/111	50/114	58/112	67/130	312	IP68 ¹	Y	Y	Y	N
Audéo M70-312T	050-0497-xx	B612PHO	C2	46/111	50/114	58/112	67/130	312	IP68 ¹	Y	Y	Y	N
Audéo M50-312T	050-0496-xx	B613PHO	C2	46/111	50/114	58/112	67/130	312	IP68 ¹	Y	Y	Y	N
Audéo M30-312T	050-0495-xx	B614PHO	C2	46/111	50/114	58/112	67/130	312	IP68 ¹	Y	Y	Y	N

* Noahlink wireless essential for programming.

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Colour chart



P1

Sand Beige



P3

Sandalwood



P4

Chestnut



P5

Champagne



P6

Silver Grey



P7

Graphite Grey



P8

Velvet Black



T7

Alpline White



H0

Beige

Phonak Audéo™ M Rechargeable

Premium (90), Advanced (70), Standard (50), Essential (30)



Audéo M-R



Audéo M-RT

Model	Order Reference Refer to colour chart	HSP Code	Category	Gain / MPO (2cc Coupler)				Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS B Compatibility
				S	M	P	UP						
Audéo M90-R*	050-0486-xx	B600PHO	C2	46/111	50/114	58/112	67/130	Rechargeable	IP68 ¹	Y	N	Y	N
Audéo M70-R*	050-0485-xx	B602PHO	C2	46/111	50/114	58/112	67/130	Rechargeable	IP68 ¹	Y	N	Y	N
Audéo M50-R*	050-0440-xx	B604PHO	C2	46/111	50/114	58/112	67/130	Rechargeable	IP68 ¹	Y	N	Y	N
Audéo M30-R*	050-0439-xx	B606PHO	C2	46/111	50/114	58/112	67/130	Rechargeable	IP68 ¹	Y	N	Y	N
Audéo M90-RT*	050-0490-xx	B607PHO	C2	46/111	50/114	58/112	67/130	Rechargeable	IP68 ¹	Y	Y	Y	N
Audéo M70-RT*	050-0489-xx	B608PHO	C2	46/111	50/114	58/112	67/130	Rechargeable	IP68 ¹	Y	Y	Y	N
Audéo M50-RT*	050-0488-xx	B609PHO	C2	46/111	50/114	58/112	67/130	Rechargeable	IP68 ¹	Y	Y	Y	N
Audéo M30-RT*	050-0487-xx	B610PHO	C2	46/111	50/114	58/112	67/130	Rechargeable	IP68 ¹	Y	Y	Y	N

* Noalink wireless essential for programming.

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Charging options

Phonak Charger BTE RIC 075-3016-10

This charger is suitable for all Paradise rechargeable BTE and RIC hearing aids. It features a holding space for large earmoulds, a USB-C charging inlet and an indicator light.



Phonak Charge and Care 075-0041-10

The Phonak Charge and Care is a charger that dries hearing aids with active ventilated air flow and treats hearing aids with UV light.



Phonak Charger Case Combi 075-0040

The Phonak Charger Case Combi is a charger, drying kit and protective hard case all in one.



Phonak Power Pack 071-0003

For charging rechargeable hearing aids on the go, the Phonak Power Pack can be easily attached to the bottom of the Phonak Charger Case Combi.



Colour chart

P1	P3	P4	P5	P6	P7	P8	T7	H0
Sand Beige	Sandalwood	Chestnut	Champagne	Silver Grey	Graphite Grey	Velvet Black	Alpine White	Beige

[Table of contents](#)

Page 12

Phonak Bolero™ M

Premium (90), Advanced (70), Standard (50), Essential (30)



Bolero M-M

Model	Order Reference Refer to colour chart	HSP Code	Category	Gain / MPO (Ear Simulator)		Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS B Compatibility
				Slimtube HE	HE11 680						
Bolero M90-M*	050-0535-xx	B635PHO	C2	58/122	63/128	312	IP68 ¹	Y	Y	Y	N
Bolero M70-M*	050-0534-xx	B636PHO	C2	58/122	63/128	312	IP68 ¹	Y	Y	Y	N
Bolero M50-M*	050-0533-xx	B637PHO	C2	58/122	63/128	312	IP68 ¹	Y	Y	Y	N
Bolero M30-M*	050-0532-xx	B638PHO	C2	58/122	63/128	312	IP68 ¹	Y	Y	Y	N

* Noahlink wireless essential for programming.

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Colour chart



P1

Sand Beige



P3

Sandalwood



P4

Chestnut



P5

Champagne



P6

Silver Grey



P7

Graphite Grey



P8

Velvet Black



T7

Alpine White



H0

Beige

Phonak Bolero™ M Rechargeable

Premium (90), Advanced (70), Standard (50), Essential (30)



Bolero M-PR

Model	Order Reference Refer to colour chart	HSP Code	Category	Gain / MPO (Ear Simulator)		Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS B Compatibility
				Slimtube HE	HE11 680						
Bolero M90-PR*	050-0677-xx	B639PHO	C2	60/127	66/130	Rechargeable	IP68 ¹	Y	N	Y	N
Bolero M70-PR*	050-0676-xx	B640PHO	C2	60/127	66/130	Rechargeable	IP68 ¹	Y	N	Y	N
Bolero M50-PR*	050-0675-xx	B641PHO	C2	60/127	66/130	Rechargeable	IP68 ¹	Y	N	Y	N
Bolero M30-PR*	050-0674-xx	B642PHO	C2	60/127	66/130	Rechargeable	IP68 ¹	Y	N	Y	N

* Noalink wireless essential for programming.

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Charging option

Phonak Charger BTE RIC 075-3016-10

This charger is suitable for all rechargeable Marvel BTE and RIC hearing aids. It features a holding space for large earmoulds, a USB-C charging inlet and an indicator light.



Colour chart

P1	P3	P4	P5	P6	P7	P8	T7	H0
Sand Beige	Sandalwood	Chestnut	Champagne	Silver Grey	Graphite Grey	Velvet Black	Alpine White	Beige

Phonak Naída™ M

Premium (90), Advanced (70), Standard (50), Essential (30)



Model	Order Reference Refer to colour chart	HSP Code	Category	Gain / MPO (Ear Simulator)		Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS B Compatibility
				SlimTube HE	HE11						
Naída M90-SP	050-0579-xx	B659PH0	C1	70/135	81/139	13	IP68 ¹	Y	Y	Y	N
Naída M70-SP	050-0578-xx	B660PH0	C1	70/135	81/139	13	IP68 ¹	Y	Y	Y	N
Naída M50-SP	050-0577-xx	B661PH0	C1	70/135	81/139	13	IP68 ¹	Y	Y	Y	N
Naída M30-SP	050-0576-xx	B662PH0	C1	70/135	81/139	13	IP68 ¹	Y	Y	Y	N

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Colour chart

P1	P3	P4	P5	P6	P7	P8	T7	H0
Sand Beige	Sandalwood	Chestnut	Champagne	Silver Grey	Graphite Grey	Velvet Black	Alphine White	Beige

Phonak Virto™ M

Premium (90), Advanced (70), Standard (50), Essential (30)




Model	Order Reference Refer to colour chart	HSP Code	Category	Gain / MPO (Ear Simulator)				Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	Microphone Technology	CROS B Compatibility
				M	P	SP	UP							
Virto M90-10 NW O	063-0490-01	C303PHO	C3	50/120	60/125	70/129	79/135	10	IP68 ¹	Y	N	N	Omni	N
Virto M70-10 NW O	063-0491-01	C305PHO	C3	50/120	60/125	70/129	79/135	10	IP68 ¹	Y	N	N	Omni	N
Virto M50-10 NW O	063-0492-01	C307PHO	C3	50/120	60/125	70/129	79/135	10	IP68 ¹	Y	N	N	Omni	N
Virto M30-10 NW O	063-0493-01	C309PHO	C3	50/120	60/125	70/129	79/135	10	IP68 ¹	Y	N	N	Omni	N
Virto M90-312 NW O	063-0494-01	C302PHO	C3	50/120	60/125	70/129	79/135	312	-	Y	N	N	Omni	N
Virto M70-312 NW O	063-0495-01	C304PHO	C3	50/120	60/125	70/129	79/135	312	-	Y	N	N	Omni	N
Virto M50-312 NW O	063-0496-01	C306PHO	C3	50/120	60/125	70/129	79/135	312	-	Y	N	N	Omni	N
Virto M30-312 NW O	063-0497-01	C308PHO	C3	50/120	60/125	70/129	79/135	312	-	Y	N	N	Omni	N
Virto M90-312	063-0498-01	I136PHO	C3	50/120	60/125	70/129	79/135	312	-	N	Y	Y	Directional	N
Virto M70-312	063-0499-01	I137PHO	C3	50/120	60/125	70/129	79/135	312	-	N	Y	Y	Directional	N
Virto M50-312	063-0500-01	I138PHO	C3	50/120	60/125	70/129	79/135	312	-	N	Y	Y	Directional	N
Virto M30-312	063-0501-01	I139PHO	C3	50/120	60/125	70/129	79/135	312	-	N	Y	Y	Directional	N

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Colour chart – Faceplate



Colour chart – Shell



Phonak Virto™ M-Titanium

Premium (90), Advanced (70)



Virto M-312

Model	Order Reference Refer to colour chart	HSP Code	Category	Gain / MPO (Ear Simulator)			Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	Microphone Technology	CROS B Compatibility
				M	P	SP							
Virto M90-Titanium	063-0488-01	C300PHO	C3	50/120	60/125	70/129	10	IP68 ¹	Y	N	N/A	Omni	N/A
Virto M70-Titanium	063-0489-01	C301PHO	C3	50/120	60/125	70/129	10	IP68 ¹	Y	N	N/A	Omni	N/A


¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Colour chart – Faceplate

Faceplate colour

- 06 Black
- 22 Cocoa
- 26 Pink

Shell colour
Titanium Grey



Wireless Communication Portfolio



MiniControl (without push button only)

Tools



Titanium FitGuide Tool
054-0847

Phonak Sky™ M

Premium (90), Advanced (70), Standard (50)




Model	Order Reference Refer to colour chart	OHS Code	OHS Category	Gain / MPO (Ear Simulator)							Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS B Compatibility
				Std xR	P xR	UP xR	SlimTube HE	HE10 680	HE11 680	HE11						
Sky M90-M	050-0648-xx	B695PHO	C2	-	-	-	58/122	63/128	-	-	312	IP68 ¹	Y	Y	Y	N
Sky M90-SP	050-0756-xx	B697PHO	C1	-	-	-	70/135	81/139	-	-	13	IP68 ¹	Y	Y	Y	N
Sky M90-PR	050-0746-xx	B696PHO	C2	-	-	-	60/127	66/130	-	-	Rechargeable	IP68 ¹	Y	Y	Y	N
Sky M70-M	050-0647-xx	B698PHO	C2	-	-	-	58/122	63/128	-	-	312	IP68 ¹	Y	Y	Y	N
Sky M70-SP	050-0755-xx	B700PHO	C1	-	-	-	70/135	81/139	-	-	13	IP68 ¹	Y	Y	Y	N
Sky M70-PR	050-0745-xx	B699PHO	C2	-	-	-	60/127	66/130	-	-	Rechargeable	IP68 ¹	Y	Y	Y	N
Sky M50-M	050-0635-xx	B701PHO	C2	-	-	-	58/122	63/128	-	-	312	IP68 ¹	Y	Y	Y	N
Sky M50-SP	050-0754-xx	B703PHO	C1	-	-	-	70/135	81/139	-	-	13	IP68 ¹	Y	Y	Y	N
Sky M50-PR	050-0744-xx	B702PHO	C1	-	-	-	60/127	66/130	-	-	Rechargeable	IP68 ¹	Y	Y	Y	N


¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Colour chart

Housing



Ear hook



Hearing Instruments: Belong Platform

Phonak Audéo™ B

Premium (90), Advanced (70), Standard (50), Essential (30)



Model	Order Reference Refer to colour chart	OHS Code	OHS Category	Gain / MPO (Ear Simulator)			Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS B Compatibility
				Std xR	P xR	UP xR						
Audéo B90-10*	050-0349-xx	B459PHO	C2	56/122	65/132	-	10	IP68 ¹	Y	N	Y	Y
Audéo B90-312	050-0353-xx	B463PHO	C2	56/122	65/132	72/136	312	IP68 ¹	Y	N	Y	Y
Audéo B90-312T	050-0357-xx	B467PHO	C2	56/122	65/132	72/136	312	IP68 ¹	Y	Y	Y	Y
Audéo B90-13	050-0361-xx	B471PHO	C2	56/122	65/132	72/136	13	IP68 ¹	Y	Y	Y	Y
Audéo B70-10*	050-0348-xx	B460PHO	C2	56/122	65/132	-	10	IP68 ¹	Y	N	Y	Y
Audéo B70-312	050-0352-xx	B464PHO	C2	56/122	65/132	72/136	312	IP68 ¹	Y	N	Y	Y
Audéo B70-312T	050-0356-xx	B468PHO	C2	56/122	65/132	72/136	312	IP68 ¹	Y	Y	Y	Y
Audéo B70-13	050-0360-xx	B472PHO	C2	56/122	65/132	72/136	13	IP68 ¹	Y	Y	Y	Y
Audéo B50-10*	050-0347-xx	B461PHO	C2	56/122	65/132	-	10	IP68 ¹	Y	N	Y	Y
Audéo B50-312	050-0351-xx	B465PHO	C2	56/122	65/132	72/136	312	IP68 ¹	Y	N	Y	Y
Audéo B50-312T	050-0355-xx	B469PHO	C2	56/122	65/132	72/136	312	IP68 ¹	Y	Y	Y	Y
Audéo B50-13	050-0359-xx	B473PHO	C2	56/122	65/132	72/136	13	IP68 ¹	Y	Y	Y	Y
Audéo B30-10*	050-0346-xx	B462PHO	C2	56/122	65/132	-	10	IP68 ¹	Y	N	Y	Y
Audéo B30-312	050-0350-xx	B466PHO	C2	56/122	65/132	72/136	312	IP68 ¹	Y	N	Y	Y
Audéo B30-312T	050-0354-xx	B470PHO	C2	56/122	65/132	72/136	312	IP68 ¹	Y	Y	Y	Y
Audéo B30-13	050-0358-xx	B474PHO	C2	56/122	65/132	72/136	13	IP68 ¹	Y	Y	Y	Y

* iCube II essential for programming

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Colour chart



Phonak Audéo™ B Rechargeable

Premium (90), Advanced (70), Standard (50)



Audéo B-R

Model	Order Reference Refer to colour chart	OHS Code	OHS Category	Gain / MPO (Ear Simulator)			Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS B Compatibility
				Std xR	P xR	UP xR						
Audéo B90-R* ¹	050-0237-xx	B475PHO	C2	56/122	65/132	72/136	Rechargeable	IP68 ²	Y	N	Y	Y
Audéo B70-R* ¹	050-0235-xx	B476PHO	C2	56/122	65/132	72/136	Rechargeable	IP68 ²	Y	N	Y	Y
Audéo B50-R* ¹	050-0234-xx	B477PHO	C2	56/122	65/132	72/136	Rechargeable	IP68 ²	Y	Y	Y	Y

* iCube II essential for programming

¹ Audéo B-R orders are currently supplied with the Phonak charging case.

² IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Charging option

Phonak Charger Case Combi 075-0040

The Phonak Charger Case Combi is a charger, drying kit and protective hard case all in one.



Phonak Power Pack 071-0003

For charging rechargeable hearing aids on the go, the Phonak Power Pack can be easily attached to the bottom of the Phonak Charger Case Combi.



Colour chart



P1

Sand Beige



P5

Champagne



P6

Silver Grey



P7

Graphite Grey



P8

Velvet Black

Phonak Bolero™ B

Premium (90), Advanced (70), Standard (50), Essential (30)



Model	Order Reference Refer to colour chart	OHS Code	OHS Category	Gain / MPO (Ear Simulator)		Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS B Compatibility
				Slimtube HE	HE10						
Bolero B90-M	050-0365-xx	B489PHO	C2	60/126	68/134	312	IP68 ¹	Y	Y	Y	Y
Bolero B90-P	050-0369-xx	B490PHO	C2	62/130	71/135	13	IP68 ¹	Y	Y	Y	Y
Bolero B90-SP	050-0373-xx	B491PHO	C2	72/132	75/137	13	IP68 ¹	Y	Y	Y	Y
Bolero B70-M	050-0364-xx	B486PHO	C2	60/126	68/134	312	IP68 ¹	Y	Y	Y	Y
Bolero B70-P	050-0368-xx	B487PHO	C2	62/130	71/135	13	IP68 ¹	Y	Y	Y	Y
Bolero B70-SP	050-0372-xx	B488PHO	C2	72/132	75/137	13	IP68 ¹	Y	Y	Y	Y
Bolero B50-M	050-0363-xx	B483PHO	C2	60/126	68/134	312	IP68 ¹	Y	Y	Y	Y
Bolero B50-P	050-0367-xx	B484PHO	C2	62/130	71/135	13	IP68 ¹	Y	Y	Y	Y
Bolero B50-SP	050-0371-xx	B485PHO	C2	72/132	75/137	13	IP68 ¹	Y	Y	Y	Y
Bolero B30-M	050-0362-xx	B480PHO	C2	60/126	68/134	312	IP68 ¹	Y	Y	Y	Y
Bolero B30-P	050-0366-xx	B481PHO	C2	62/130	71/135	13	IP68 ¹	Y	Y	Y	Y
Bolero B30-SP	050-0370-xx	B482PHO	C2	72/132	75/137	13	IP68 ¹	Y	Y	Y	Y

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Colour chart



P1

Sand Beige



P5

Champagne



P6

Silver Grey



P7

Graphite Grey



P8

Velvet Black

Phonak Bolero™ B Rechargeable

Premium (90), Advanced (70), Standard (50)



Bolero B-PR

Model	Order Reference Refer to colour chart	OHS Code	OHS Category	Gain / MPO (Ear Simulator)		Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS B Compatibility
				Slimtube HE	HE10						
Bolero B90-PR* ¹	050-0376-xx	B494PHO	C2	62/130	71/135	Rechargeable	IP68 ²	Y	Y	Y	N
Bolero B70-PR* ¹	050-0375-xx	B493PHO	C2	62/130	71/135	Rechargeable	IP68 ²	Y	Y	Y	N
Bolero B50-PR* ¹	050-0374-xx	B492PHO	C2	62/130	71/135	Rechargeable	IP68 ²	Y	Y	Y	N

* iCube II essential for programming

¹ Audéo B-R orders are currently supplied with the Phonak charging case.

² IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Charging option

Phonak Charger Case Combi 075-0040

The Phonak Charger Case Combi is a charger, drying kit and protective hard case all in one.



Phonak Power Pack 071-0003

For charging rechargeable hearing aids on the go, the Phonak Power Pack can be easily attached to the bottom of the Phonak Charger Case Combi.



Colour chart



P1

Sand Beige



P5

Champagne



P6

Silver Grey



P7

Graphite Grey




P8

Velvet Black

Phonak CROS B

Compatible with Belong™ Platform



Model	Order Reference Refer to colour chart	HSP Code	Battery Size	IP Rating	Compatibility
CROS B-312 Custom	063-0246-01	N0187CR	312	IP68 ¹	Audéo B and Bolero B*
CROS B-13 Custom	063-0247-01	N0188CR	13	IP68 ¹	Audéo B and Bolero B*
CROS B-312	050-0274-xx	N0179CR	312	IP68 ¹	Audéo B and Bolero B*
CROS B-13	050-0275-xx	N0180CR	13	IP68 ¹	Audéo B and Bolero B*

* Not compatible with Audeo B-R devices
¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

CROS B performance

CROS B also offers the Speech in Loud Noise program with StereoZoom. It is available in combination with B90 and B70 performance level hearing aids. For the lower performance levels, CROS B provides Real Ear Sound and omni directional microphone modes. On the better ear, CROS B clients fully benefit from all Belong features like AutoSense OS. It provides seamless transitions and optimum settings for different listening situations.

Hearing aid performance level	Available microphone modes for CROS B			
	Automatic StereoZoom	Manual StereoZoom	Real Ear Sound	Omni directional
Premium (B90)	•	•	•	•
Advanced (B70)		•	•	•
Standard (B50)			•	•
Essential (B30)				•

Colour chart

Housing

P1 Sand Beige	P5 Champagne	P6 Silver Grey
P7 Graphite Grey	P8 Velvet Black	

Shell

14 Tan	22 Cocoa	26 Pink	28 Brown

Faceplate

21 Clear Transparent	14 Tan	22 Cocoa	26 Pink	28 Brown

Earpiece options

CROS demo hook	- Hook to hold CROS on the ear during fitting and demonstration - For demo purposes only - Included with all Phonak CROS devices
CROS SlimTube	- Standard SlimTube for CROS devices - Available in left/right and four different lengths (0, 1, 2, 3) - Domes need to be attached. Available as open, closed and power dome in three sizes (S, M, L)
CROS Hook	- Earhook to attach a custom earmould with standard tubing to CROS -13 models - Only available for the CROS -13 models - Included with any CROS -13 device or can be ordered separately (order number: 048-0261)
CROS Tip	- Custom made earpiece for CROS devices - Available in left/right and four different lengths (0, 1, 2, 3)

Phonak CROS B Rechargeable

Compatible with Belong™ Platform



Model	Order Reference Refer to colour chart	HSP Code	Battery Size	IP Rating	Compatibility
CROS B-R*	050-0601-xx	N0201CR	Rechargeable	IP68 ¹	Audéo B-R

* iCube II essential for programming
¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Colour chart



P1
Sand Beige



P5
Champagne



P6
Silver Grey



P7
Graphite Grey



P8
Velvet Black

Phonak Naída™ B

Premium (90), Advanced (70), Standard (50), Essential (30)



Model	Order Reference Refer to colour chart	OHS Code	OHS Category	Gain / MPO (Ear Simulator)		Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS B Compatibility
				HE11 680	HE11						
Naída B90-SP	050-0398-xx	B548PHO	C1	77/139	83/143	13	IP68 ¹	Y	Y	Y	Y
Naída B90-UP	050-0404-xx	B569PHO	C1	81/140	85/145	675	IP68 ¹	Y	Y	Y	Y
Naída B70-SP	050-0397-xx	B567PHO	C1	77/139	83/143	13	IP68 ¹	Y	Y	Y	Y
Naída B70-UP	050-0403-xx	B570PHO	C1	81/140	85/145	675	IP68 ¹	Y	Y	Y	Y
Naída B50-SP	050-0319-xx	B568PHO	C1	77/139	83/143	13	IP68 ¹	Y	Y	Y	Y
Naída B50-UP	050-0402-xx	B571PHO	C1	81/140	85/145	675	IP68 ¹	Y	Y	Y	Y
Naída B30-SP	050-0318-xx	B572PHO	C1	77/139	83/143	13	IP68 ¹	Y	Y	Y	Y
Naída B30-UP	050-0399-xx	B573PHO	C1	81/140	85/145	675	IP68 ¹	Y	Y	Y	Y

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Colour chart



P1

Sand Beige



P5

Champagne



P6

Silver Grey



P7

Graphite Grey



P8

Velvet Black

Phonak Naída™ B Rechargeable

Premium (90), Advanced (70), Standard (50)



Naída B-R

Model	Order Reference Refer to colour chart	OHS Code	OHS Category	Gain / MPO (Ear Simulator)			Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS B Compatibility
				Std xR	P xR	UP xR						
Naída B90-R RIC ^{*1}	050-0437-xx	B545PHO	C2	56/112	65/132	72/136	Rechargeable	IP68 ²	Y	N	Y	N
Naída B70-R RIC ^{*1}	050-0436-xx	B546PHO	C2	56/112	65/132	72/136	Rechargeable	IP68 ²	Y	N	Y	N
Naída B50-R RIC ^{*1}	050-0429-xx	B547PHO	C2	56/112	65/132	72/136	Rechargeable	IP68 ²	Y	N	Y	N

* iCube II essential for programming

¹ Naída B-R orders are currently supplied with the Phonak charging case.

² IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Charging option

Phonak Charger Case Combi 075-0040

The Phonak Charger Case Combi is a charger, drying kit and protective hard case all in one.



Phonak Power Pack 071-0003

For charging rechargeable hearing aids on the go, the Phonak Power Pack can be easily attached to the bottom of the Phonak Charger Case Combi.



Colour chart



P1

Sand Beige



P5

Champagne



P6

Silver Grey



P7

Graphite Grey



P8

Velvet Black

Phonak Naída™ Link



Model	Order Reference Refer to colour chart	OHS Code	Category	Gain / MPO (Ear Simulator)						Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS II Compatibility
				St dxR	P xR	UP xR	HE7	HE7 680	Power Slimtube						
Naída Link-RIC	050-0262-xx	B479PHO	C2	121/56	132/65	137/72	-	-	-	13	IP67 ¹	N	Y	Y	N
Naída Link-UP	050-0261-xx	B478PHO	C1	-	-	-	144/85	139/80	144/85	675	IP67 ¹	N	Y	Y	N

¹ IP67 indicates that the hearing aid is water resistant and dust tight. It survived immersion in 1 metre of water for 30 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Phonak Sky™ B

Premium (90), Advanced (70), Standard (50)



Model	Order Reference Refer to colour chart	OHS Code	OHS Category	Gain / MPO (Ear Simulator)							Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS B Compatibility
				St dxR	P xR	UP xR	SlimTube HE	HE10 680	HE11 680	HE11						
Sky B90-M	050-0465-xx	B549PHO	C2	-	-	-	126/60	134/68	-	-	312	IP68 ¹	Y	Y	Y	Y
Sky B90-P	050-0469-xx	B552PHO	C2	-	-	-	130/62	135/71	-	-	13	IP68 ¹	Y	Y	Y	Y
Sky B90-SP	050-0473-xx	B558PHO	C1	-	-	-	-	-	139/77	143/83	13	IP68 ¹	Y	Y	Y	Y
Sky B90-UP	050-0477-xx	B561PHO	C1	-	-	-	-	-	140/81	145/85	675	IP68 ¹	Y	Y	Y	Y
Sky B90-RIC	050-0481-xx	B564PHO	C2	122/56	132/65	136/72	-	-	-	-	13	IP68 ¹	Y	Y	Y	Y
Sky B70-M	050-0464-xx	B550PHO	C2	-	-	-	126/60	134/68	-	-	312	IP68 ¹	Y	Y	Y	Y
Sky B70-P	050-0468-xx	B553PHO	C2	-	-	-	130/62	135/71	-	-	13	IP68 ¹	Y	Y	Y	Y
Sky B70-SP	050-0472-xx	B559PHO	C1	-	-	-	-	-	139/77	143/83	13	IP68 ¹	Y	Y	Y	Y
Sky B70-UP	050-0476-xx	B562PHO	C1	-	-	-	-	-	140/81	145/85	675	IP68 ¹	Y	Y	Y	Y
Sky B70-RIC	050-0480-xx	B565PHO	C2	122/56	132/65	136/72	-	-	-	-	13	IP68 ¹	Y	Y	Y	Y
Sky B50-M	050-0463-xx	B551PHO	C2	-	-	-	126/60	134/68	-	-	312	IP68 ¹	Y	Y	Y	Y
Sky B50-P	050-0467-xx	B554PHO	C2	-	-	-	130/62	135/71	-	-	13	IP68 ¹	Y	Y	Y	Y
Sky B50-SP	050-0471-xx	B560PHO	C1	-	-	-	-	-	139/77	143/83	13	IP68 ¹	Y	Y	Y	Y
Sky B50-UP	050-0475-xx	B563PHO	C1	-	-	-	-	-	140/81	145/85	675	IP68 ¹	Y	Y	Y	Y
Sky B50-RIC	050-0479-xx	B566PHO	C2	122/56	132/65	136/72	-	-	-	-	13	IP68 ¹	Y	Y	Y	Y

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Colour chart

Housing



Ear hook

R1 Green	R2 Yellow	R3 Pink
R4 Orange	13 Transparent	Q6 Blue
Q7 Purple		

Phonak Sky™ B Rechargeable

Premium (90), Advanced (70), Standard (50)



Model	Order Reference Refer to colour chart	OHS Code	OHS Category	Gain / MPO (Ear Simulator)		Battery Size	IP Rating	Tinnitus Balance	Telecoil	Wireless	CROS B Compatibility
				SlimTube HE	HE10 680						
Sky B90-PR* ¹	050-0484-xx	B555PHO	C2	130/62	135/71	Rechargeable	IP68 ²	Y	Y	Y	N
Sky B70-PR* ¹	050-0483-xx	B556PHO	C2	130/62	135/71	Rechargeable	IP68 ²	Y	Y	Y	N
Sky B50-PR* ¹	050-0482-xx	B557PHO	C2	130/62	135/71	Rechargeable	IP68 ²	Y	Y	Y	N

* iCube II essential for programming

¹ Naida B-R orders are currently supplied with the Phonak charging case.

² IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Charging option

Phonak Charger Case Combi 075-0040

The Phonak Charger Case Combi is a charger, drying kit and protective hard case all in one.



Phonak Power Pack 071-0003

For charging rechargeable hearing aids on the go, the Phonak Power Pack can be easily attached to the bottom of the Phonak Charger Case Combi.



Phonak Virto™ B

Premium (90), Advanced (70), Standard (50), Essential (30)



Model	Order Reference Refer to colour chart	OHS Code	Category	Gain / MPO (Ear Simulator)				Battery Size	Telecoil	Wireless	Microphone Technology	CROS B Compatibility
				M	P	SP	UP					
Virto B90-10 NW O	063-0206-01	C262PHO	C3	50/120	60/125	70/129	79/135	10	Y	N	Omni	N
Virto B90-312 NW O	063-0287-01	C274PHO	C3	50/120	60/125	70/129	79/135	312	Y	N	Omni	N
Virto B90-10 O	063-0211-01	C266PHO	C3	50/120	60/125	70/129	79/135	10	Y	Y	Omni	N
Virto B90-10	063-0216-01	C270PHO	C3	50/120	60/125	70/129	79/135	10	Y	Y	Directional	N
Virto B90-312	063-0159-01	C278PHO	C3	50/120	60/125	70/129	79/135	312	Y	Y	Directional	Y
Virto B90-13	063-0236-01	I126PHO	C3	50/120	60/125	70/129	79/135	13	Y	Y	Directional	Y
Virto B70-10 NW O	063-0207-01	C263PHO	C3	50/120	60/125	70/129	79/135	10	Y	N	Omni	N
Virto B70-312 NW O	063-0288-01	C275PHO	C3	50/120	60/125	70/129	79/135	312	Y	N	Omni	N
Virto B70-10 O	063-0212-01	C267PHO	C3	50/120	60/125	70/129	79/135	10	Y	Y	Omni	N
Virto B70-10	063-0217-01	C271PHO	C3	50/120	60/125	70/129	79/135	10	Y	Y	Directional	N
Virto B70-312	063-0224-01	C279PHO	C3	50/120	60/125	70/129	79/135	312	Y	Y	Directional	Y
Virto B70-13	063-0237-01	I127PHO	C3	50/120	60/125	70/129	79/135	13	Y	Y	Directional	Y
Virto B50-10 NW O	063-0208-01	C264PHO	C3	50/120	60/125	70/129	79/135	10	Y	N	Omni	N
Virto B50-312 NW O	063-0482-01	C276PHO	C3	50/120	60/125	70/129	79/135	312	Y	N	Omni	N
Virto B50-10 O	063-0213-01	C268PHO	C3	50/120	60/125	70/129	79/135	10	Y	Y	Omni	N
Virto B50-10	063-0218-01	C272PHO	C3	50/120	60/125	70/129	79/135	10	Y	Y	Directional	N
Virto B50-312	063-0234-01	C280PHO	C3	50/120	60/125	70/129	79/135	312	Y	Y	Directional	Y
Virto B50-13	063-0238-01	I128PHO	C3	50/120	60/125	70/129	79/135	13	Y	Y	Directional	Y
Virto B30-10 NW O	063-0209-01	C265PHO	C3	50/120	60/125	70/129	79/135	10	Y	N	Omni	N
Virto B30-312 NW O	063-0483-01	C277PHO	C3	50/120	60/125	70/129	79/135	312	Y	N	Omni	N
Virto B30-10 O	063-0214-01	C269PHO	C3	50/120	60/125	70/129	79/135	10	Y	Y	Omni	N
Virto B30-10	063-0219-01	C273PHO	C3	50/120	60/125	70/129	79/135	10	Y	Y	Directional	N
Virto B30-312	063-0235-01	C281PHO	C3	50/120	60/125	70/129	79/135	312	Y	Y	Directional	Y
Virto B30-13	063-0239-01	I129PHO	C3	50/120	60/125	70/129	79/135	13	Y	Y	Directional	Y

Hearing Instruments: Free-to-Client

Phonak Vitus

Free-to-Client Devices



Model	Order Reference Refer to colour chart	OHS Code	Category	Gain / MPO (Ear Simulator)				Battery Size	IP Rating	Wireless	CROS B Compatibility
				Slim Tube	HE7	HE7 680	Power Slimtube				
Vitus BTE-micro	050-0573-xx	B544PH0	C2	60/126	-	60/131	-	312	IP68 ¹	N	N
Vitus BTE-P	050-0240-xx	B543PH0	C2	62/170	-	68/138	-	13	IP68 ¹	N	N
Vitus BTE-UP	050-0269-xx	B542PH0	C1	81/140	-	73/139	-	675	IP68 ¹	N	N

¹ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 metre of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.

Colour chart



Hearing Protection

Phonak Serenity Choice™

Generic hearing protection for occasional use

These protective earplugs come in three standard sizes, with different filters for each application.

A wide range of applications: music, motorsport, work, shooting & hunting, comfort, sleep.



Serenity Choice	Order Reference
Serenity Choice Music	108-3085-006010
Serenity Choice Sleep	108-3086-006002
Serenity Choice Motorsport	108-3084-006002
Serenity Choice Work	108-3089-006002
Serenity Choice Comfort	108-3087-006012
Serenity Choice Hunting & Shooting	108-3088-006006
Serenity Choice Starter Kit	108-3100

Serenity Choice generic hearing protection has a minimum order quantity of 10.

Not returnable for credit or exchange

Serenity Filter Kits	Filter Type	Order Reference
Serenity Choice Motorsport + (S) or Serenity Choice Work + (A) (S)	Filter Phonak K120 Kit	004-4408-003
Serenity Choice Motorsport + (S) or Serenity Choice Work + (A) (S)	Filter Phonak K125 Kit	004-4409-003
Serenity Choice Work + (A) (S)	Filter Phonak KI30 Kit	004-4410-003
Serenity Choice Shooting + Hunting + (A) (S)	Filter Phonak KIM9 Kit	004-4414-003
Serenity Choice Comfort + (S)	Filter Phonak KM10 Kit	004-4415-003
Serenity Choice Music + (A) (S)	Filter Phonak KM15 Kit	004-4416-003
Serenity Choice Fly + (A) (S)	Filter Phonak KM16 Kit	004-4551-003
Serenity Choice Music + (A) (S)	Filter Phonak KM20 Kit	004-4417-003
Serenity Choice Music + (A) (S)	Filter Phonak KM25 Kit	004-4418-003
Serenity Choice Swim + (S)	Filter Phonak KR5 Kit	004-4420-003
NA	Pressor with tag clear kit	005-3232-3

One kit contains a package of 10 units



Serenity Choice Starter Kit

Phonak Serenity Choice™ Plus

Custom made hearing protection for regular use

Custom earplugs with different filters for each application; individually sculpted for a perfect fit and maximum hearing protection.

Available in two materials (silicone and acrylic) and two size options (Full Concha and Mini).

A wide range of applications: music, motorsport, work, shooting & hunting, fly, comfort, sleep, swim.




Hearing protection	Order Reference	Material Type
Serenity Choice Plus Motorsport + (S)	108-3106	Silicone
Serenity Choice Plus Fly + (A)	108-3044	Acrylic
Serenity Choice Plus Fly + (S)	108-3066	Silicone
Serenity Choice Plus Music + (A)	108-3041	Acrylic
Serenity Choice Plus Music + (S)	108-3067	Silicone
Serenity Choice Plus Work + (A)	108-3037	Acrylic
Serenity Choice Plus Work + (S)	108-3107	Silicone
Hearing comfort		
Serenity Choice Plus Sleep + (S)	108-3045	Silicone
Serenity Choice Plus Swim + (S)	108-3038	Silicone
Serenity Choice Plus Comfort + (S)	108-3110	Silicone
Serenity Choice Plus Shooting + Hunting + (A)	108-3108	Acrylic
Serenity Choice Plus Shooting + Hunting + (S)	108-3109	Silicone



Not returnable for credit or exchange

Colour chart

Soft Silicone



Hard Acrylic



Mini



Full Concha

Programming Tools

Programming Tools

Programmes	Order Reference
Phonak Target™ Fitting Software	Contact Phonak Customer Service
NOAHLink Wireless	058-0139
NOAH 4 Software - Users License Pack	058-5151
Hi-Pro 2 Hardware	058-0110
LisN-S	052-3414

Wireless Fitting Devices	Order Reference
iCube II	076-0026-P809
iCube II charger	075-3000-20
iCube II Spare Neckloop	076-0029



iCube II and Accessory Set



NOAHLink Wireless

Programming Cables	Order Reference	Product Compatibility
Plug-in Cable - 10A Long Left (2m)	058-0098	Tao 10, Virto B10, Virto Titanium
Plug-in Cable - 10A Long Right (2m)	058-0099	Tao 10, Virto B10, Virto Titanium
Plug-in Cable - 10A Short Left (0.5m)	058-1071	Tao 10, Virto B10, Virto Titanium
Plug-in Cable - 10A Short Right (0.5m)	058-1070	Tao 10, Virto B10, Virto Titanium
Plug-In Cable 312 - Left	058-0030	Custom products using 312 battery
Plug-In Cable 312 - Right	058-0029	Custom products using 312 battery
Plug-In Cable 13 - Left	058-0043	Custom products using 13 battery
Plug-In Cable 13 - Right	058-0042	Custom products using 13 battery
10A Programming Cable Kit (Long 2m)	058-1079	Virto B10, Virto Titanium
Programming Cable 10A - Left (0.5m)	058-1071	Virto / Nano
Programming Cable 10A - Right (0.5m)	058-1070	Virto / Nano
Silicon Ear Guide for 10A Programming Cables	004-2311	Nano
4 pin BTE Cable - Left (Long - 1m)	058-0256	BTE, RIC
4 pin BTE Cable - Right (Long - 1m)	058-0257	BTE, RIC
4 pin BTE Cable, Short CS 44 (0.5m)	058-0255	BTE, RIC
Mini USB Cable	076-0853	
Micro USB Cable	075-3003	
Phonak programming cable sleeves	004-8220-001	

Wireless Accessories

Phonak Wireless Communication Portfolio

Streamers	Order Reference	Compatible with
ComPilot II	076-0023-P520	Venture, Belong

Transmitter	Order Reference	Compatible with
RemoteMic II	076-0039-P509	ComPilot Air II, ComPilot II, ComPilot
TVLink II	076-0040-P520	ComPilot Air II, ComPilot II, ComPilot
PartnerMic™	076-4001-P510	Paradise and Marvel platform devices

Remote Control	Order Reference	Compatible with
PilotOne II	076-0019-P5	Venture, Belong
MiniControl	054-1006	Virto B-Titanium, Virto M-Titanium, Virto B-10 NWO, Virto M-10 NWO
ComPilot II	076-0023-P520	Venture, Belong
RemoteControl	076-0065-P5	Paradise and Marvel platform devices

Television	Order Reference	Compatible with
TV Connector	076-3002-06	Paradise, Marvel, Audeo B-Direct



ComPilot II



Remote Mic II



TVLink II



PilotOne II



MiniControl



TV Connector



RemoteControl



PartnerMic™

* GST will be added to invoice

† Non-returnable item

Phonak Wireless Communication Portfolio

Accessories

Accessories	Order Reference	Compatible with
Replacement Antenna Neckloop (Short)	076-0012 - Black; 076-0034 - Grey	ComPilot, ComPilot II
Replacement Antenna Neckloop (Long)	076-0013 - Black; 076-0035 - Grey	ComPilot, ComPilot II
Silicone Protection Cover	005-0369	ComPilot, ComPilot II
Cinch Adaptor	052-0371	Audio Cable
Audio Cable 1.5m (3.5mm jack male to cinch)	052-0423	RCA connections on older televisions
Audio Input Cable 3.5mm (1m)	052-0359	All accessories
USB Charger/Plug (select cable separately)	075-3000-20	Roger Pen, Roger EasyPen, Roger Clip-On Mic, Roger Touchscreen Mic
Micro USB Cable for Charger/Plug	075-3003	USB Charger/Plug
Mini USB Cable for Charger/Plug	075-0853	USB Charger/Plug
TosLink Cable	057-0011	TV Connector, Roger Pen, Roger EasyPen, Roger Clip-On Mic, Roger Touchscreen Mic
External Lapel Microphone for ComPilot	057-0115	ComPilot
USB-A to USB-C cable 1m flexible	043-3093	Charger BTE RIC Case, Mini Charger Case RIC, Partner Mic

Roger™ for Marvel and Paradise

Roger Microphones	Order Reference (refer to colour chart)
Roger On iN (microphone only)	056-3024-xx010
Roger Select iN (analogue docking station)	056-3002-xx
Roger Table Mic II iN	056-4001-P5
Roger Table Mic II iN (x2)	056-4002-P5

Roger Installation	Order Reference (refer to colour chart)
Roger Installer	055-3033
Roger X (03)	052-3215-X03P5
Roger X (02)	052-3113-X02P5

Roger iN: How does it work?

Roger provides exceptional hearing performance over distance and in noise using its tried and tested transmitter and receiver system.

Roger iN transmitters are provided with two on board receiver installations which are easily transferred to a pair of Paradise or Marvel devices.

All other Roger transmitters can still be used with Paradise or Marvel devices but virtual receivers need to be transferred from a Roger X via Roger Installer.

Colour chart

Roger On iN

P5	P7	P7	P5	W3
Champagne	Graphite Grey	Graphite Grey	Champagne	Pearl White



Roger On iN Roger Table Mic II iN Roger Select iN Roger Installer

Roger™ Portfolio

Roger Receivers (03) Compatible with Roger Clip-on Mic, Roger Select, Roger On, Roger TableMic II

Design Integrated Receivers for Phonak Hearing Aids	Order Reference (refer to colour chart)	Compatible with
Roger 19 (03)	052-3292-D03xx	Naida B-UP, Naida V-UP, Sky V-UP
Roger 18 (03)	052-3283-D03xx	Naida B-SP, Audéo B-13, Naida V-SP, Naida V-RIC, Bolero V-P, Bolero V-SP, Audéo V-13, Bolero B, Sky V-P/SP

Design Integrated Receivers for Implants	Order Reference (refer to colour chart)	Compatible with
Roger 17 (03)	052-3236-D03xx	Advanced Bionics: Naida CI
Roger 14 (03)	052-3223-D03xx	Cochlear: Nucleus 5, Nucleus 6

[View colour chart](#)

Design Integrated Receivers for Implants	Order Reference (refer to colour chart)	Compatible with
Roger X (03)	052-3215-X03P5	Must be used with an AudioShoe (AS) or ComPilot/ComPilot II Streamer for Phonak Hearing aids and AB Implants
AS19	054-0732-xx	Naida B-UP, Naida V-UP, Sky V-UP
AS18	054-0244-xx	Naida B-SP, Audéo B-13, Naida V-SP, Naida V-RIC, Bolero V-P, Bolero V-SP, Audéo V-13, Bolero B-13, Sky V-P/SP
AS16	054-0119-xx	Bolero Q-P, Bolero Q-SP
AS13	054-0536-xx	Ambra microP/SP, Solana microP/SP, Cassia microP/SP, Dalia microP/SP
Roger NeckLoop (03)	056-4005-P5110	

Adaptors	Order Reference	Compatible with
Roger 20 Adapter	055-3008	Roger 20 Receiver
Roger 21 Adapter	055-3020	Roger 21 Receiver



Roger 19



Roger 18



Roger 17



Roger 14



Roger X

Roger™ Portfolio

Roger Microphones

Roger Microphones	Order Reference (refer to colour chart)
Roger On (microphone only)	056-3023-xx010
Roger Clip-On Mic	052-3294
Roger Table Mic II	052-4118-000P5
Roger Select (analogue docking station)	052-3471-000-xx

Bundles	Order Reference (refer to colour chart)
Roger TableMic II Bundle	052-4119-000P5

[View colour chart](#)

Accessories	Order Reference	Compatible with
Roger Clip-On Mic Docking/Charging Station	055-3023	Roger Clip-On Mic
Roger Select Docking/Charging Station	055-3025	Roger Select
Neck Loop Lanyard	052-3198	
Neck Loop Lanyard	052-3234	Roger Clip-On Mic
Storage Case	017-3275	
Pouch	017-3278	Roger Clip-On Mic
Remote Control	055-4075-T4	Roger Table Mic II
Soft Case - Black	017-4004	Roger Table Mic II
USB Charger/Plug (select cable separately)	075-3000-20	Roger Clip-On Mic, Roger Touchscreen Mic
Micro USB Cable for Charger/Plug	075-3003	USB Charger/Plug
Mini USB Cable for Charger/Plug	076-0853	USB Charger/Plug
Roger Select Clip	055-3030	Roger Select
Lanyard Roger Select	055-3029	Roger Select
Roger Table Mic II Remote Control	055-4075-T4	Roger Table Mic II



Roger On



Roger
Clip-On Mic



Roger
Table Mic II



Roger Table Mic II
Bundle



Roger Select

Roger™ Portfolio for Education

Roger Receivers (02) and Microphones

Design Integrated Receivers for Phonak Hearing Aids	Order Reference (refer to colour chart)	Compatible with
Roger 19 (02)	052-3291-D02xx	Sky V-UP, Naída V-UP
Roger 18 (02)	052-3281-D02xx	Sky V-SP, Sky V-P, Sky V-RIC, Naída V-SP, Naída V-RIC, Bolero V-P, Bolero V-SP, Audéo V-13

Design Integrated Receivers for Implants	Order Reference (refer to colour chart)	Compatible with
Roger 21 (02)	052-3469-D02xx	Sonnet
Roger 20 (02)	052-3317-D02xx	Nucleus 7
Roger 17 (02)	052-3235-D02xx	Advanced Bionics: Naída CI
Roger 14 (02)	052-3222-D02xx	Cochlear: Nucleus 5, Nucleus 6

[View colour chart](#)

Universal Roger Receivers	Order Reference (refer to colour chart)	Compatible with
Roger X (02)	052-3113-X02P5	Used with DAI port
Roger NeckLoop (02)	056-4004-P5210	
Roger Focus II	056-3006-xx	
Roger Focus II-312	056-3005-xx	

Roger Transmitters	Order Reference (refer to colour chart)	Compatible with
Roger Touchscreen Mic	052-3297	Roger (02) Receivers, DSF Speakers and Roger for Education Microphones
Roger Pass-around	052-3261	
Roger Multimedia Hub	052-3299	

Other	Order Reference	Compatible with
Roger Charging Rack	055-3010	Roger for Education Transmitters
Roger Repeater	052-3476	All Roger for Education devices
Roger BaseStation	052-3342	All Roger for Education devices



Roger™ Dynamic Soundfield

Receivers and Transmitters

Receivers	Order Reference
Roger DigiMaster 5000 loudspeaker	052-3195
Roger DigiMaster 7000 loudspeaker	052-3197
Roger WallPilot	052-3115

Receivers Accessories	Order Reference
Wall Mount Kit	052-3428
Floor Stand Complete	052-3164
Tripod Stand for Digimaster Speaker	052-3444
Tool Kit	026-3008
MLx Audio Checker Set	052-3177
Adapter MLx BAHA	052-3002
ML 14i Checker	052-3434

Transmitters	Order Reference
Roger BaseStation	052-3342

Transmitters Accessories	Order Reference
Roger Charging Rack	055-3010
Carry Case for Roger Touchscreen Mic	017-3272
Lanyard for Roger Touchscreen Mic	055-3005
Protection cover for Roger Touchscreen Mic	005-3184-006
Boom mic 3.5mm for Roger Touchscreen Mic	055-3009
MiniUSB Programming Cable	043-3027
Audio Cable 1 metre	052-0359
Audio Cable 3 metres	052-0340
Housing Set for Roger WallPilot	052-3242
Cinch adapter for Audio Cable	052-0371
Belt clip for Roger Touchscreen Mic (to be used with Boom mic)	004-4031



Roger™ Dynamic Soundfield

Power Cables, Bundles and Kits

Power Cables	Order Reference
Power Supply for DigiMaster loudspeaker	052-3141
Power Cable 2pin 1.8m for Digimaster loudspeaker	052-3149
USB Charger/Plug (select cable separately)	075-3000-20
Micro USB Cable for Charger/Plug	075-3003
Mini USB Cable for Charger/Plug	076-0853

Bundle Packages	Order Reference
Educational Classroom Kit	052-3464
DSF 5000 with Roger Touchscreen Starter Kit - Floor Stand	052-3465
DSF 5000 with Roger Touchscreen Starter Kit - Wall Mount	052-3466
DSF 7000 with Touchscreen Starter Kit - Wall Mount	052-3467
DSF 7000 Starter Kit with Touchscreen - Floor Stand	052-3468

Educational Classroom Kit
1 x Roger Touchscreen Mic 1 x Roger Multimedia Hub 2 x Roger Pass-around 1 x Roger Charging Rack

DSF 5000 Roger Touchscreen Kit
1 x DigiMaster 5000 LoudSpeaker 1 x Roger Touchscreen 1 x Floor Stand (base + tube) or Wall Mount or Tripod Stand 1 x Tool Kit 1 x Power Adaptor and Cable 1 x Charger 1 x User Guides

DSF 7000 Roger Touchscreen Kit
2 x DigiMaster 7000 LoudSpeaker 1 x Roger Touchscreen 2 x Floor Stand (base + tube) or Wall Mount or Tripod Stand 1 x Tool Kit 2 x Power Adaptor and Cable 1 x Charger 1 x User Guides

Phonak Colour Chart

Roger Receiver (02/03) and Audio Shoe Colours

	P1 Sand Beige	P2 Amber Beige	P3 Sandalwood	P4 Chestnut	P5 Champagne	P6 Silver Grey	P7 Graphite Grey	P8 Velvet Black	P9 Ruby	Q1 Petrol	Q2 Electric Green	Q3 Caribbean Pirate	T3 Precious Pink	O1 Beige	L8 White	L9 Black	M1 Beige/Sand	M2 Charcoal	T1 Mocha/Brown	T2 Smoke/Grey	XN/T7 Alpine White	XP Princess Pink	M6 Lava Red	M7 Blue Ocean	M8 Majesty Purple	V1 Anthracite	V2 Beige	V3 Black	V4 Ebony	V5 Nordic Grey	V6 White
Roger 17	●			●		●		●	●	●											●	●									
Roger 18	●				●	●	●	●				●	●										●	●	●						
Roger 19	●				●	●	●	●				●	●										●	●	●						
Roger 20															●	●	●		●	●											
Roger 21																										●	●	●	●	●	●

Roger On Colours

P5	P7	P7	P5	W3
Champagne	Graphite Grey	Graphite Grey	Champagne	Pearl White

Roger Focus Receiver Colours

P1	P2	P3	P4	P6	P8	P9	Q1	Q2	Q3
Sand Beige	Amber Beige	Sandalwood	Chestnut	Silver Grey	Velvet Black	Ruby	Petrol	Electric Green	Caribbean Pirate
Q4	Q5	Q6	Q7	T3	O1	13			
Dragon Orange	Vanity Pink	Blue Lagoon Transparent	Royal Purple Transparent	Precious Pink	Beige	Pure Transparent			

Roger Focus II & Roger Focus II-312* Receiver Colours

M6	M7	M8	Q2	Q3	T3	T7	P6	P8	H0
Lava Red	Ocean Blue	Majesty Purple	Electric Green	Caribbean Pirate	Precious Pink	Alpine White	Silver Grey*	Velvet Black*	Beige

Sound Delivery Systems

SDS 4.0 Receivers & SlimTubes

Paradise and Marvel RIC and BTE models



Standard 4.0 Receivers		Compatible with	Pack Qty
054-0777 0L	054-0778 0R	Paradise and Marvel	1
054-0779 1L	054-0780 1R	Paradise and Marvel	1
054-0781 2L	054-0782 2R	Paradise and Marvel	1
054-0783 3L	054-0784 3R	Paradise and Marvel	1

Moderate 4.0 Receivers		Compatible with	Pack Qty
054-0801 0L	054-0802 0R	Paradise and Marvel	1
054-0803 1L	054-0804 1R	Paradise and Marvel	1
054-0805 2L	054-0806 2R	Paradise and Marvel	1
054-0807 3L	054-0808 3R	Paradise and Marvel	1

Power 4.0 Receivers		Compatible with	Pack Qty
054-0812 0L	054-0813 0R	Paradise and Marvel	1
054-0814 1L	054-0815 1R	Paradise and Marvel	1
054-0816 2L	054-0817 2R	Paradise and Marvel	1
054-0818 3L	054-0819 3R	Paradise and Marvel	1

Receiver Tools 4.0	Order Reference	Compatible with	Pack Qty
Receiver Handling Tool for Slim Tips 4.0	026-1528	Paradise and Marvel	1
Measurement tool for receiver length	026-0697	Paradise and Marvel	1
Pin Removal Tool	002-0404-1	Paradise and Marvel	1
CeruShield™ Disk	098-0445	Paradise and Marvel	8


Receiver Retention 4.0	Order Reference	Compatible with	Pack Qty
Retention for S Receiver	004-8596-001	Paradise and Marvel	10
Retention for M Receiver	004-8597-001	Paradise and Marvel	10
Retention for P Receiver	004-8623-001	Paradise and Marvel	10

SlimTubes	Order Reference	Compatible with	Pack Qty
SlimTube 4.0 00L	054-0884	Bolero M	5
SlimTube 4.0 00R	054-0885	Bolero M	5
SlimTube 4.0 0L	054-0794	Bolero M	5
SlimTube 4.0 0R	054-0795	Bolero M	5
SlimTube 4.0 1L	054-0796	Bolero M	5
SlimTube 4.0 1R	054-0798	Bolero M	5
SlimTube 4.0 2L	054-0799	Bolero M	5
SlimTube 4.0 2R	054-0823	Bolero M	5
SlimTube 4.0 3L	054-0824	Bolero M	5
SlimTube 4.0 3R	054-0825	Bolero M	5

Fit'nGo Kits	Order Reference	Compatible with
RIC Fit'nGo Kit 4.0	054-0879	RIC Paradise and Marvel
SlimTube Fit'nGo Kit 4.0	054-0894	BTE Marvel

SDS 4.0 Domes & Custom Coupling

Paradise and Marvel RIC and BTE models



Domes	Order Reference	Compatible with	Pack Qty
Cap Dome 4.0 (S)	054-0788	Paradise and Marvel	10
Open Dome 4.0 (S)	054-0785	Paradise and Marvel	10
Open Dome 4.0 (M)	054-0786	Paradise and Marvel	10
Open Dome 4.0 (L)	054-0787	Paradise and Marvel	10
Vented Dome 4.0 (S)	054-0809	Paradise and Marvel	10
Vented Dome 4.0 (M)	054-0810	Paradise and Marvel	10
Vented Dome 4.0 (L)	054-0811	Paradise and Marvel	10
Power Dome 4.0 (S)	054-0820	Paradise and Marvel	10
Power Dome 4.0 (M)	054-0821	Paradise and Marvel	10
Power Dome 4.0 (L)	054-0822	Paradise and Marvel	10

Custom Coupling	Order Reference	Compatible with
SlimTip 4.0 Acrylic	062-0006-01	Paradise and Marvel
SlimTip 4.0 Silicone	062-0007-01	Paradise and Marvel
cShell 4.0 Acrylic	062-0008-01	Paradise and Marvel
cShell 4.0 Titanium	062-0009-01	Paradise and Marvel*
CROS Tip 4.0 Acrylic	062-0016-**	Paradise

* Chargeable with aid order – not complimentary.

Not returnable for credit or exchange

Receivers, Domes & Coupling

Belong, Venture and older models only



Standard Receivers	Order Reference	Compatible with	Pack Qty
xReceiver Unit xS 0-R	054-0103	Spice, Quest, Venture, Belong	1
xReceiver Unit xS 1-R	054-0105	Spice, Quest, Venture, Belong	1
xReceiver Unit xS 2-R	054-0107	Spice, Quest, Venture, Belong	1
xReceiver Unit xS 3-R	054-0109	Spice, Quest, Venture, Belong	1
xReceiver Unit xS 0-L	054-0102	Spice, Quest, Venture, Belong	1
xReceiver Unit xS 1-L	054-0104	Spice, Quest, Venture, Belong	1
xReceiver Unit xS 2-L	054-0106	Spice, Quest, Venture, Belong	1
xReceiver Unit xS 3-L	054-0108	Spice, Quest, Venture, Belong	1

Power Receivers	Order Reference	Compatible with	Pack Qty
xReceiver Unit xS 0-R	054-0236	Spice, Quest, Venture, Belong	1
xReceiver Unit xP 1-R	054-0238	Spice, Quest, Venture, Belong	1
xReceiver Unit xP 2-R	054-0240	Spice, Quest, Venture, Belong	1
xReceiver Unit xP 3-R	054-0242	Spice, Quest, Venture, Belong	1
xReceiver Unit xP 0-L	054-0235	Spice, Quest, Venture, Belong	1
xReceiver Unit xP 1-L	054-0237	Spice, Quest, Venture, Belong	1
xReceiver Unit xP 2-L	054-0239	Spice, Quest, Venture, Belong	1
xReceiver Unit xP 3-L	054-0241	Spice, Quest, Venture, Belong	1

Domes	Order Reference	Compatible with	Pack Qty
Open Domes (S)	054-1987	Spice, Quest, Venture, Belong	10
Open Domes (M)	054-1988	Spice, Quest, Venture, Belong	10
Open Domes (L)	054-1989	Spice, Quest, Venture, Belong	10
Closed Domes (S)	054-1990	Spice, Quest, Venture, Belong	10
Closed Domes (M)	054-1991	Spice, Quest, Venture, Belong	10
Closed Domes (L)	054-1992	Spice, Quest, Venture, Belong	10
Power Domes (S)	054-1993	Spice, Quest, Venture, Belong	10
Power Domes (M)	054-1994	Spice, Quest, Venture, Belong	10
Power Domes (L)	054-1995	Spice, Quest, Venture, Belong	10

Fit'nGo Kits	Order Reference	Compatible with
SlimTube HE Fit'nGo Kit	054-0314	BTE
SlimTube Fit'nGo Kit	054-0316	BTE
RIC Fit'nGo Kit	054-0087	RIC Quest, Venture, Belong

Custom Coupling	Order Reference	Compatible with
cShell	062-0014-02	Receiver
Slim Tip	062-0012-02	Slim Tube, Receiver
CROS Tip	062-0025-01	CROS, CROS II, CROS B

Receivers, Domes & Coupling

Belong, Venture and older models only

Earmould Material	Order Reference
Hard Acrylic BTE Standard Earmould	-
Silicone 40 BTE Standard Earmould	-
Silicone 70 BTE Standard Earmould	-
Earmould	062-0001-01

SlimTubes



SlimTubes (Blister Packs)	Order Reference	Compatible with	Pack Qty
SlimTube 0B-R	054-0324	CORE, Baseo Q	5
SlimTube 1B-R	054-0325	CORE, Baseo Q	5
SlimTube 2B-R	054-0326	CORE, Baseo Q	5
SlimTube 3B-R	054-0327	CORE, Baseo Q	5
SlimTube 0B-L	054-0328	CORE, Baseo Q	5
SlimTube 1B-L	054-0329	CORE, Baseo Q	5
SlimTube 2B-L	054-0330	CORE, Baseo Q	5
SlimTube 3B-L	054-0331	CORE, Baseo Q	5

SlimTubes HE (Blister Packs)	Order Reference	Compatible with	Pack Qty
SlimTube HE 0-R	054-0539	Spice, Quest, Venture, Belong	5
SlimTube HE 1-R	054-0540	Spice, Quest, Venture, Belong	5
SlimTube HE 2-R	054-0541	Spice, Quest, Venture, Belong	5
SlimTube HE 3-R	054-0542	Spice, Quest, Venture, Belong	5
SlimTube HE 0-L	054-0543	Spice, Quest, Venture, Belong	5
SlimTube HE 1-L	054-0544	Spice, Quest, Venture, Belong	5
SlimTube HE 2-L	054-0545	Spice, Quest, Venture, Belong	5
SlimTube HE 3-L	054-0546	Spice, Quest, Venture, Belong	5

Power SlimTube II (Blister Packs)	Order Reference	Compatible with	Pack Qty
Power SlimTube II 00-R	054-0755	Naída V/B-SP/UP, Sky V/B-SP/UP, Vitus UP	2
Power SlimTube II 0-R	054-0756	Naída V/B-SP/UP, Sky V/B-SP/UP, Vitus UP	2
Power SlimTube II 1-R	054-0757	Naída V/B-SP/UP, Sky V/B-SP/UP, Vitus UP	2
Power SlimTube II 2-R	054-0758	Naída V/B-SP/UP, Sky V/B-SP/UP, Vitus UP	2
Power SlimTube II 3-R	054-0759	Naída V/B-SP/UP, Sky V/B-SP/UP, Vitus UP	2
Power SlimTube II 00-L	054-0760	Naída V/B-SP/UP, Sky V/B-SP/UP, Vitus UP	2
Power SlimTube II 0-L	054-0761	Naída V/B-SP/UP, Sky V/B-SP/UP, Vitus UP	2
Power SlimTube II 1-L	054-0762	Naída V/B-SP/UP, Sky V/B-SP/UP, Vitus UP	2
Power SlimTube II 2-L	054-0763	Naída V/B-SP/UP, Sky V/B-SP/UP, Vitus UP	2
Power SlimTube II 3-L	054-0764	Naída V/B-SP/UP, Sky V/B-SP/UP, Vitus UP	2

CROS SlimTube (Blister Packs)	Order Reference	Compatible with	Pack Qty
CROS SlimTube 0-R	054-0686	CROS, CROS II, CROS B	1
CROS SlimTube 1-R	054-0687	CROS, CROS II, CROS B	1
CROS SlimTube 2-R	054-0688	CROS, CROS II, CROS B	1
CROS SlimTube 3-R	054-0689	CROS, CROS II, CROS B	1
CROS SlimTube 0-L	054-0690	CROS, CROS II, CROS B	1
CROS SlimTube 1-L	054-0691	CROS, CROS II, CROS B	1
CROS SlimTube 2-L	054-0692	CROS, CROS II, CROS B	1
CROS SlimTube 3-L	054-0693	CROS, CROS II, CROS B	1

Tone Hooks



Hook HE11 680/
Hook HE11



Hook HE10 680/
Hook HE10



Mini Hook HE10
680/Mini Hook HE10



Hook HE9 680



Mini Hook HE9 680/
Mini Hook HE9



Hook HE7 680/
Hook HE7



Mini Hook HE7 680/
Mini Hook HE7



CROS II Hook

Tone Hooks	Order Reference	Product Compatibility	Pack Qty
Hook HE11 680	042-0098-131	Naida B, Naida V, Sky V, Bolero M, Naida M, Sky M	10
Hook HE 11	004-1667-131	Naida B, Naida V, Sky V, Bolero M, Naida M, Sky M	10
Hook HE10 680	004-8028-001	Ambra, Solana, Cassia, Bolero Q, Bolero V & B	10
Hook HE 10	004-1471-001	Ambra, Solana, Cassia, Bolero Q, Bolero V & B	10
Mini Hook HE10 680	004-8041-001	Ambra, Solana, Cassia, Bolero Q, Bolero V & B	10
Mini Hook HE10	004-1472-001	Ambra, Solana, Cassia, Bolero Q, Bolero V & B	10
Hook HE7 680	004-1207-001	Ambra, Solana, Cassia, Bolero Q, Bolero V & B	10
Hook HE7	004-1865-001	Exélia Art, Versáta Art, Certéna Art, Eleva, eXtra, Una, Savia Art, Naida UP, Naida SP, Baseo Q	10
Mini Hook HE7 680	004-2775-001	Exelia Art, Versáta Art, Certéna Art Eleva, eXtra, Una, Savia Art, Naida UP, Naida SP, Baseo Q	10
Mini Hook HE7	004-1230-001	Naida UP, Naida SP, Baseo Q	10
CROS II Hook Set	048-0261	CROS II-13, CROS B-13 <i>Includes: Pin removal tool, spare pin, instruction guide</i>	1

Fit'nGo Kits	Order Reference	Compatible with
Hook Fit'nGo kit	054-0750	Sky V

- **HE 680 (Horn Effect 860 Ohm)** is an optimal tone hook and gives overall smoothing of the resonant peaks with minimal gain loss. Less suppression at 1000Hz is offered compared to the PD 1000.
- **HE (Horn Effect)** tone hooks provide extended high frequency delivery and are generally used as the standard hook on high gain models.
- **Mini Tone Hooks** offer the same acoustic characteristics as the standard range but are specially designed for users with small ears.
- Ear hooks with no numbers are undamped and the numbers refer to the amount of damping.
- CROS II Hook is used to attach a custom earmould with standard tubing to a CROS II-13 model only.

Cleaning and Wax Management

Cleaning Kit	Order Reference	Product Compatibility
D-Dry+	098-0511-00014	All models
C&C Line Cleaning Kit 1	098-0527	BTE
C&C Line Cleaning Kit 2	098-0528	micro, ITE

Cleaning Kit Components	Order Reference	Pack Qty
C&C Cleansing Tablets	098-0530	1
C&C Cleansing Spray 8ml	098-0529	1
C&C Cleansing Tissue Refill	098-0526	25
C&C Drying Capsules Refill	098-0014	6
Air Puffer	098-0018	1
Drying/Cleansing Beaker	098-0194	1

Wax Management	Order Reference	Pack Qty
CeruShield™ Disk	098-0445	8^
Filter HF4 blue (Dispenser)	098-0374	15
Filter HF4 red (Dispenser)	098-0373	15
Filter HF3 blue (Dispenser)	098-0293	15
Filter HF3 red (Dispenser)	098-0292	15
Cerustop	098-0282	8 ⁺
Wax Brush	098-0064	1

^10 packets per box
 +100 packets per box



D-Dry+

C&C Line Cleaning Kit

Not returnable for credit or exchange

Microphone Protection

Microphone Protection	Order Reference	Product Compatibility	Pack QTY
Microphone Protection	004-1234-131	MicroBTE	1
Microphone Protection	042-0038-011	M13, Spice H20, Naída S CRT	1
Microphone Protection	004-1868-011	Baseo Q	1
Microphone Port Cover - Omni	044-0161-26	Custom	1
MicGuard	054-0082-13	Nano	8
Wind/weather Protection	004-1281-001		1

Batteries

Private Label^ Mercury-Free PowerOne	Order Reference	Battery Size
PowerOne	070-0364	10
PowerOne	070-0365	312
PowerOne	070-0366	13
PowerOne	070-0367	675


Mercury-Free PowerOne	Order Reference	Battery Size
PowerOne	070-0293	10
PowerOne	070-0295	312
PowerOne	070-0294	13
PowerOne	070-0296	675

Battery Accessories	Order Reference
Battery tester & battery holder keyring	098-0347

Batteries can only be ordered in multiples of 300 cells ie: per carton.



PowerOne p10
Mercury-Free



PowerOne p312
Mercury-Free



PowerOne p13
Mercury-Free



PowerOne p675
Mercury-Free

^For Private Label artwork, please contact your local Phonak Representative.

Not returnable for credit or exchange

Consumables

Ear Impression Materials	Order Reference
Otoform Kc - Impression material with hardner (750g)	095-0003
Otoform Ak - A+B 2 component impression material (400g in each tub)	095-0005

Earmould Tubing	Order Reference	Qty
Sound tube, 13 medium tube	004-8227-21	25
Sound tube, 13 thick wall tube	004-8228-21	25
Sound tube, 13 dry wall tube	004-8229-21	10
Sound tube, 13 dry thick wall tube	004-8230-21	10
TRS-13 Clear Thick Tube w/Tube Lock	021-5277	10
TRS-13 Clear Thick Dri Tube w/Tube Lock	004-6376	10

Ear Impression Essentials	Order Reference	Qty
Heine Otoscope	054-0286	1
Heine Ear Light	026-1150	1
Heine Otoscope Bulbs 2000 (S)	026-5646	
Heine Otoscope Bulbs 3000 (L)	026-5647	
Heine Otoscope Tips (2.5mm)	019-0223	50
Heine Otoscope Tips (4mm)	019-0224	50
Impression Syringe	026-5654	
Tubing Expanders	026-5009	
Phonak Impression boxes (for shipping)	017-5255	50

Ear Impression Essentials	Order Reference	Qty
EasyView Otoblocks^ S	054-0769	20
EasyView Otoblocks^ M	054-0768	20
EasyView Otoblocks^ L	054-0767	20
Kot-N-Blox - Blue (S)	095-5025	50
Kot-N-Blox - White (M)	095-5026	50
Kot-N-Blox - Red (L)	095-5027	50

Tools

Description	Order Reference	Min. Qty
Tweezers for soldering work	026-0032	
Tweezers Ideal 21	026-0040	
Screwdriver for Housing Screw (1.3mm blade)	026-0010	
Screwdriver (2.00mm)	026-0013	
Phonak Screwdriver	004-0235	
Holder for Scalpel	026-0025	
Side Cutter	026-0015	
Single Edge Blades	026-5179	100
Surgical Scissors	026-5099	
ITE Battery Door Extraction Tool	026-5264	
Unlocking tool (size 13 RIC removal tool)	003-0251-1	
Titanium FitGuide Tool	054-0847	
Small Vent Cleaning Tool	026-1415	
Cleaning Tool	098-0064	
Power Supply AU/NZ	075-0020-14	

Description	Order Reference	Min. Qty
Pin Removal Tool	002-0404-1	
Stethoclip	021-0306	
ITE attachment for Stethoclip	021-0447	
Cleaning Pins (for Slim Et micro Tubes)	054-0344	
Drill bit for Removal Pin	026-5648	
Tubing Expanders	026-5009	
SlimTip Removal Tool	054-0402	
Flexfit 3 Transparent	005-0474-021	
Flexfit Tube 3 Transparent	005-0475-021	
Lubricant Earmould Auragel	098-0379	10
Magnet (for use with phone)	054-5117	

Policies and purchase conditions

Policies and purchase conditions

Price

The prices in this product catalogue are subject to change at the discretion of Phonak with or without notice to you. Hearing instruments are exempt from GST in Australia. The GST applicable for all other products are listed in the product catalogue. If GST is payable, it will be added to your Tax Invoice. A Tax Invoice will be emailed to the designated account email address.

Orders

Any order not accompanied by a current order form may experience delays in processing.

All orders must be emailed to orders@phonak.com.au, posted to **12 Inglewood Place, Norwest NSW 2153** or faxed to **1800 212 515**, using a current Phonak order form.

Product order forms can be requested via the above email address, by calling 1800 212 313 or by visiting the [Phonak Australia website](https://www.phonak.com.au).

Phonak takes all reasonable endeavours to meet your requested time frame for the delivery of orders once the order form is processed. The delivery time may be impacted if you do not complete the mandatory fields on the order form including but not limited to the 'Required By' section. Phonak will endeavour to deliver your orders to you by the times stipulated below. Note, these times exclude transit time.

eStore

To order online via the eStore visit Phonak Pro, select Store in the top banner and sign in with your PhonakPro ID. Available via this link PhonakPro.com

Internal turn-around times*			
Despatch (business days)			
Order	Standard	VIP	Comments for customers
BTE/RIC New orders (with no earmoulds)	2	n/a	Standard turnaround times for orders received prior to 2pm EST on a business day^
Earmoulds and Custom Devices New BTE/RIC orders with earmoulds New orders Remakes	10* 10* 10*	2 2 2	Turnaround times for orders received: Standard: Prior to 10am EST on a business day VIP: Repairs Prior to 12pm EST on a business day
Repairs Custom Devices BTE/RIC	5 5	2 2	
Other HI consumables Accessories Batteries Customer labelled batteries	2 2 2 on request	n/a n/a n/a n/a	Standard stock items Standard stock items Standard stock items

*General observations. If these turnaround times cannot be met, you will be contacted by Phonak Customer Care.

Phonak offers a 48hr express service for all repairs and new earmould and custom device orders (IIC orders and non-standard colours may require longer processing time). This service will carry a fee and must be requested at time of order (\$55 incl GST per order).

^If 'Required By' date is stipulated on your order form, your order will be dispatched to meet this date and as such, may not be released within the times stipulated in the table above.

Warranty

Hearing Instrument Warranty

Hearing instruments purchased from Phonak are supplied with an international 3 year manufactures' warranty, covering manufacturing and material defects, commencing on the date of the original Tax Invoice. This warranty does not cover damage caused by improper handling, wax or moisture damage. Repairs made by unauthorised parties are not covered and will void the warranty.

Receivers purchased with, and at the time of hearing instrument order come with a 3 year manufacturers' warranty. Receivers purchased singularly, outside of hearing instrument warranty come with a 1 year manufacturers' warranty.

An additional one year warranty (up to a maximum of 4 year warranty) is available to purchase within 90 days of the hearing instrument purchase for the charges indicated below.

Free-To-Client: \$85 per device (plus GST)

Top Up: \$160 per device (plus GST)

Earmoulds and any customised ear tips and / or shells on custom hearing instruments are supplied with a 90 day warranty, such that any ear mould / shell remakes with the same specifications, required within the first 3 months of invoice, will be made at no charge. Any remake that involves an alteration of the original order specification will incur a charge and is not covered by the warranty. All earmoulds are not returnable for credit.

Accessory Warranty

Remote Controls and Accessories are supplied with a 12 month limited manufacturers' warranty, covering manufacturing and material defects. Additional warranties are not available for accessories. Proof of purchase must be provided with claim.

Rechargeable cases are supplied with a 3 year manufacturer's warranty.

International Warranty

A 12 month international warranty is provided on all serialised Hearing Instruments, Audiometric equipment and Roger products which arise when the user is out of the original country of purchase. Warranty is valid on presentation of warranty card or invoice issued by authorised clinic.

Loss and damage

Phonak hearing instruments can be supplied with a one-time only loss and damage coverage for 12 months from the date of original invoice for fees indicated below. No substitutions accepted. A 'Loss or Damage Claim Form – Statutory Declaration' must be completed to claim a replacement and submitted to Phonak. One claim can be made per invoiced instrument that is damaged or lost within the first 12 months of the date of the original Tax Invoice. Subsequent loss or damage of the replacement instrument does not attract any coverage. The warranty duration of the originally supplied instrument is maintained (i.e. no new warranty is applied to the replacement device). Instruments replaced under loss and damage are not eligible for credit.

Loss and damage replacement fee

550

Return for credit

Free-to-Client hearing instruments are accepted for exchange only within the 60 day exchange period from the date of invoice.

Free-to-client hearing instruments are not accepted for credit.

All Top-Up hearing instruments are accepted for return for credit within 90 days from the date of original invoice.

Batteries, Phonak D-Dry kits or consumable stock items are not returnable for credit or exchange.

All products exchanged or returned must be returned to Phonak in their original, undamaged packaging and accompanied with a completed Credit Return Request Form. Damaged or unusable packaging may incur a fee.

Payment terms

Payment terms are 30 days from the date of statement. Acceptable forms of payment include Direct Deposit (preferred), Visa and MasterCard.

Delivery

All orders, including accessories, under \$100 will incur a \$10 (plus GST) shipping fee.

Bulk orders may take longer to arrive and, at our discretion, may incur a delivery charge which will be at your expense.

An authorised signatory must sign for the products upon delivery. The products will only be delivered to your registered business address.

Service and repair charges

Cost and inclusions for repairs:

- (i) Clean & Check: Includes full clean, inside and out and a full check of the device performance performance. A clean and check charge is only applied if no fault is found with the devices.
- (ii) RIC & BTE repairs: Includes replacement of full electronic parts. An additional 12 months warranty is applied post chargeable repair.
- (iii) ITE Repairs: Includes repair or electronics. An additional 12 months warranty is applied post chargeable repair.
- (iv) ITE re-shell: Excludes repair of electronics. An additional 3 months remake warranty is applied post chargeable repair.

Hearing instrument repairs

Within the hearing instrument warranty period
(unless the damage is caused by improper handling, wax or moisture) **No charge**

Any unauthorised repairs made to the products will automatically void the warranty. A change in vent size, style, material, colour and/or customer damage to a remake order will be deemed a new order and applicable charges will apply. The GST payable will be added to the Tax Invoice.

Outside of the warranty period but less than 5 years old the following charges apply:

Clean and Check	\$85 (incl GST)
RIC/BTE/ITE repairs	\$185 (incl GST)
ITE re-shell	\$110 (incl GST)
ITE re-shell and repair	\$295 (incl GST)

Remakes

Custom shells – under a 90 day warranty	No charge
Custom shells – outside of warranty (classed as new order)	\$95 (incl GST)
SlimTips – outside of warranty (classed as new order)	\$35 (incl GST)
Earmould re-tubing	\$15 (incl GST)

Digital Wireless Accessories (DWA) repairs

Accessories (including remote controls)	\$185 incl GST
Battery replacements	\$49.50 incl GST

Please note:

- (a) All hearing instruments must be deemed serviceable or they will be considered Damaged Beyond Repair;
- (b) There is no guarantee that obsolete hearing instruments over 5 years old can be repaired, re-shelled or remade;
- (c) Repairs on custom instruments do not include custom shell; and
- (d) Replaced parts under any chargeable repair will receive 1 year warranty on those parts.

All repairs outside of warranty will automatically proceed without quoting and charged the fixed service price as selected on the Phonak Service form.

Trial instruments

All trial instruments include a 3 year warranty on the instrument and a complimentary clean and check available at anytime.

Please note:

- Every trial instrument receives a complimentary coupling option (receiver or tube/dome)
- No additional warranty available for purchase
- Trial orders are dispatched within standard BTE shipping lead times
- No exchange or return for credit
- No loss and damage policy
- No repair service outside of warranty period
- Firmware update will be completed on trial instruments if necessary

Sale of instruments

Phonak is not open to the public. Any requests for purchase or repair from a consumer will be redirected to an authorised health care professional. The Customer commits that all end user Hearing Instrument purchasers can obtain an in-person fitting with pre-audiological and post-audiological services provided by your company, unless such in-person fitting is not required per the product's instructions for use.

Privacy

By creating an account with Phonak you accept that Phonak will store your information. We recognise that your privacy and security are of the utmost importance and Phonak is committed to securing your personal information safely. If you have any questions/concerns about our [Privacy Policy](#) please contact us at privacy@sonova.com.

Sonova Australia Pty Ltd
12 Inglewood Place,
Norwest, NSW 2153
Australia

P: 1800 809 321
E: orders@phonak.com.au
www.phonakpro.com.au