

Phonak

Sales aid

Roger Focus II

Attractive and effective solution for children and teens with unilateral hearing loss or a functional hearing loss

Who is it for?

Children or teens with:

- Unilateral hearing loss (UHL)
- Functional hearing loss resulting from:
 - Auditory Processing Disorder (APD)
 - Autism Spectrum Disorder (ASD)

What is it?

- Roger ear-level receiver
- Improves signal-to-noise ratio (SNR) and therefore the wearer's ability to understand speech at a distance and in noisy environments
- Worn in the ear(s) with normal hearing
- Fit unilaterally for UHL
- Fit bilaterally or unilaterally for APD or ASD

Unique benefits

- Only stand-alone ear-level receiver that uses Phonak Roger proprietary protocol, making it different from other FM and remote microphone systems available
 - Adaptive gain model that maintains the SNR even when the noise levels increase
 - MultiTalker network
- Used in combination with any Roger microphone
- One click and one-time pairing, along with lockable volume control, makes it simple for the child, parents and teachers to use.

