

Phonak TargetMatch support.

Phonak TargetMatch is an automated guidance system for real ear measurements (REMs). It guides you through a step-by-step workflow of calibration, probe tube placement, unaided measurements, aided measurements, and automatic target matching.

TargetMatch requires 3 components:

- Phonak Target version 9.0 or higher
- Noah software 4.4 version or higher
- REM system

You must have one of the following real ear systems to perform REMs:

- Auditdata Primus Pro software version 4.1.0 or higher
- Interacoustics Affinity Suite software version 2.19 or higher
- Natus Otosuite software version 4.88 or higher
- Signia Unity 3 software version 5.9.0 or higher

With TargetMatch, Phonak supports you in performing verification in an efficient manner.

For additional support in using TargetMatch, please contact your local Phonak representative:

Name: _____

Email: _____

Number: _____

For further guidance and support in using your REM equipment, please contact your local REM representative:

Name: _____

Email: _____

Number: _____