

Phonak CROS L Order Form

Phonak CROS L provides a CROS and BICROS solution for people with an unaidable hearing loss in one ear. CROS L is worn on the unaidable ear to pick up sounds and wirelessly transmits them to a compatible Phonak hearing aid worn on the ear with better hearing. CROS L has no acoustical output.

CROS L is compatible with Audéo L-R, Audéo L-RT, Naída L-PR and Sky L-PR hearing aids.

Step 1: Order details

Account Number: _____
 Hearing Clinic: _____
 Delivery Address: _____
 Contact Name: _____ Phone Number: _____
 Email Address: _____
 Date Ordered: _____ Date Required: _____

Step 2: Patient information

Name: _____

Step 3: VIP service

VIP 24 hour service (\$57.75 incl GST)

Step 4: Choose hearing instrument

(see following page for colour options)

Model	Part No.	HSP	Colour	Qty
CROS L-R	050-1013-xx	N0304CR		

Step 5: Choose your Acoustic Coupling Option

CROS Slim Tube 4.0 with Domes


CROS Slim Tube 4.0 is available in left/right and four different lengths (0, 1, 2, 3) with standard dome available in three sizes (S, M, L)



L R Select options in section a and c only below

CROS Slim Tube 5.0 with Slim Tip*

CROS Slim Tube 5.0 is available in left/right and four different lengths (0, 1, 2, 3), with a custom made earpiece



L R Select options in section c and d only below
 ^ Impression required

CROS Tip 4.0*

Fully custom made earpiece for CROS, available in four different cable lengths (0, 1, 2, 3)



L R Select options in section b and d only below
 ^ Impression required

Step 6: Choose from the following options

(a) Domes 4.0

Dome style	Dome size	
	Left	Right
Cap	Left	Right
Open	Left	Right
Vented	Left	Right
Power	Left	Right



(b) CROS Tip Cables – must be paired with a CROS Tip below

0-L	0-R
1-L	1-R
2-L	2-R
3-L	3-R

(c) CROS Slim Tube 5.0

054-1180	CROS Tube 5 0L	054-1184	CROS Tube 5 0R
054-1181	CROS Tube 5 1L	054-1185	CROS Tube 5 1R
054-1182	CROS Tube 5 2L	054-1186	CROS Tube 5 2R
054-1183	CROS Tube 5 3L	054-1187	CROS Tube 5 3R

(d) Slim Tip 4.0 or CROS Tip

	Shell colour	Clear (21)	Rose Tint (7)	Beige/Pink (26)	Tan (14)	Vent style	Cavity	Regular	Shell attribute	Hollow	Solid
SlimTip Acrylic (062-0006-01)											-
SlimTip Silicone (062-0007-01)											-
CROS Tip 4.0 Hard Acrylic (062-0016-**) - must be paired with a CROS Tip cable											-

Generic options (all earpieces)

Removal line Transparent (21) Removal line extension + 5mm
 Locks Canal Skeleton
 Vent size AOV
 L: ___mm R: ___mm

S = Standard, will default if no other selection ticked

Phonak CROS L Order Form

Phonak CROS L provides a CROS and BICROS solution for people with an unaidable hearing loss in one ear. CROS L is worn on the unaidable ear to pick up sounds and wirelessly transmits them to a compatible Phonak hearing aid worn on the ear with better hearing. CROS L has no acoustical output.

CROS L is compatible with Audéo L-R, Audéo L-RT, Naída L-PR and Sky L-PR hearing aids.

Step 7: Choose your wireless Lumity hearing instrument for the better ear

Audéo L-RT (NOTE: Please also complete the Phonak Audéo L-RIC order form to specify correct inclusions).

Audéo L-R (NOTE: Please also complete the Phonak Audéo L-RIC order form to specify correct inclusions).

Naída L-PR (NOTE: Please also complete the Phonak Naída L order form to specify correct inclusions).

Sky L-PR (NOTE: Please also complete the Phonak Sky L order form to specify correct inclusions).

Step 8: Choose your accessories

Complimentary charging accessory for CROS L-R
Phonak Charger Ease
(075-0080-10)

Wireless Communication Accessories
Please use the Phonak Audéo L-RIC, Naída L or Sky L order forms to specify accessories.

Roger Accessories
Refer to Phonak Roger Order Form

Step 9: Check List

Phonak Audéo L RIC, Naída L or Sky L order form accompanies this order

Roger order form accompanies this order (if applicable)

Additional information

HSP codes

Model	HSP codes	Category
CROS L-R	N0304CR	NS

Colour chart



Models



CROS L-R

Charging accessory



Phonak Charger Ease
075-0080-10