

















Phonak platform comparison.

	Features	Infinio		Lumity	Paradise	Marvel
		Sphere	R			
Phonak SmartSpeech™ Technology	 Adaptive Phonak Digital (APD) Phonak proprietary fitting formula	3.0	3.0	2.0	2.0	1.0
	 AutoSense OS Automatic operating system blending together the optimal settings for each unique sound environment	6.0	6.0	5.0	4.0	3.0
	 Spheric Speech Clarity DNN-based sound processing that separates speech from noise in any direction	Yes	-	-	-	-
	 StereoZoom Seamlessly adjusts for more speech focus from the front	2.0	2.0	2.0	1.0	1.0
	 SpeechSensor Improves hearing from the sides and back in noisy environments	90/70	90	90	-	-
	 Motion Sensor Hearing** Automatically adjusts to optimise hearing while walking	90/70	90/70/50	90/70/50	90/70/50	-
	 Dynamic Noise Cancellation Automatically reduces background noise around you	90/70	90/70	90/70	90	-
	 Speech Enhancer Boosts soft speech and distant speech in quiet	90	90	90	90	-

Phonak platform comparison.

Features		Infinio		Lumity	Paradise	Marvel
		Sphere	R			
Connectivity	 Active Bluetooth® connections Number of Bluetooth devices that can be used simultaneously	2	2	2	2	1
	 Bluetooth pairings Maximum number of paired Bluetooth devices	8	8	8	8	2
	 Auracast™ ready ERA chip	Yes	Yes	-	-	-
	 myPhonak Adjust your hearing aids to suit your needs	Yes	Yes	Yes	Yes	Yes
	 Health data tracking** Track activity levels, steps, and heart rate*	Yes	Yes	Yes	Yes	-
	 Tap Control** Pause/resume streaming, answer/end calls, and access voice assistant	90/70	90/70	90/70	90/70	-
	 RogerDirect™ Direct streaming from Roger microphones	Yes	Yes	Yes	Yes	Yes
	 CROS The Phonak solution for unaidable hearing loss in one ear	-	Yes	Yes	Yes	-

* Heart rate tracking is a feature of Audeo Fit only, requires registration

** Rechargeable models only

Some features shown are not available in Virto NW 0 & Virto Titanium hearing aids.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Sonova AG is under license.

The Auracast™ word mark and logos are trademarks owned by Bluetooth SIG, Inc. Any use of such marks by Sonova AG is under license.