

Phonak CROS L provides a CROS and BICROS solution for people with an unaidable hearing loss in one ear. CROS L is worn on the unaidable ear to pick up sounds and wirelessly transmits them to a compatible Phonak hearing aid worn on the ear with better hearing. CROS L has no acoustical output.

CROS L is compatible with Audéo L-R, Audéo L-RT, Naída L-PR and Sky L-PR hearing aids.

Step 1: Order details

Account Number:

Hearing Clinic:

Delivery Address:

Contact Name:Phone Number:

Email Address:

Date Ordered:Date Required:

Step 2: Patient information

Name:

Step 3: VIP service

VIP 24 hour service (\$57.75 incl GST)

Step 4: Choose hearing instrument (see following page for colour options)

Model	Part No.	HSP	Colour	Qty
CROS L-R	050-1013-xx	N0304CR		

Step 5: Choose your Acoustic Coupling Option

CROS Slim Tube 5.0 with Domes


CROS Slim Tube 5.0 is available in left/right and four different lengths (0, 1, 2, 3) with standard dome available in three sizes (S, M, L)



L R Select options in section a and c only below

CROS Slim Tube 5.0 with Slim Tip*

CROS Slim Tube 5.0 is available in left/right and four different lengths (0, 1, 2, 3), with a custom made earpiece



L R Select options in section c and d only below
^ Impression required

CROS Tip 4.0*

Fully custom made earpiece for CROS, available in four different cable lengths (0, 1, 2, 3)



L R Select options in section b and d only below
^ Impression required

Step 6: Choose from the following options

(a) Domes 4.0

Dome style			Dome size		
Cap	Left	Right			
Open	Left	Right	Small	Medium	Large
Vented	Left	Right	Small	Medium	Large
Power	Left	Right	Small	Medium	Large



(b) CROS Tip Cables – must be paired with a CROS Tip below

0-L 0-R
1-L 1-R
2-L 2-R
3-L 3-R

(c) CROS Slim Tube 5.0

054-1180	CROS Tube 5 0L	054-1184	CROS Tube 5 0R
054-1181	CROS Tube 5 1L	054-1185	CROS Tube 5 1R
054-1182	CROS Tube 5 2L	054-1186	CROS Tube 5 2R
054-1183	CROS Tube 5 3L	054-1187	CROS Tube 5 3R

(d) Slim Tip 4.0 or CROS Tip 4.0

	Shell colour	Clear (21)	Rose Tint (11)	Beige/Pink (26)	Tan (14)	Vent style	Cavity	Regular	Shell attribute	Hollow	Solid
SlimTip Acrylic (062-0006-01)											-
SlimTip Silicone (062-0007-01)				-	-		-			-	
CROS Tip 4.0 Hard Acrylic (062-0016-**) - must be paired with a CROS Tip cable								-			-

S = Standard, will default if no other selection ticked

Generic options (all earpieces)

Removal line Transparent (21) Removal line extension + 5mm

Locks Canal Skeleton

Vent size AOV

L:___mm R:___mm

Step 7: Choose your wireless Lumity hearing instrument for the better ear

Audéo L-RT (NOTE: Please also complete the Phonak Audéo L-RIC order form to specify correct inclusions).
Audéo L-R (NOTE: Please also complete the Phonak Audéo L-RIC order form to specify correct inclusions).
Naída L-PR (NOTE: Please also complete the Phonak Naída L order form to specify correct inclusions).
Sky L-PR (NOTE: Please also complete the Phonak Sky L order form to specify correct inclusions).

Step 8: Choose your accessories

Complimentary charging accessory for CROS L-R Phonak Charger Ease (075-0080-10)	Wireless Communication Accessories Please use the Phonak Audéo L-RIC or Naída L order forms to specify accessories.	Roger Accessories Refer to Phonak Roger Order Form
--	---	--

Step 9: Check List

Phonak Audéo L RIC or Naída L order form accompanies this order
Roger order form accompanies this order (if applicable)

Additional information

WorkSafe code

Model	HSP codes	Category
CROS L-R	HAPH282	NS

Colour chart

										
P1 Sand Beige	P3 Sandalwood	P4 Chestnut	P5 Champagne	P6 Silver Gray	P7 Graphite Gray	P8 Velvet Black	H0 Beige	Q2 Electric Green	Q3 Caribbean Pirate	T3 Precious Pink

Models



Charging accessory

