

Phonak hardware compatibility.

Our latest hearing aids

	Audéo Infinio RIC hearing aids	CROS Infinio Transmitter	Virto Infinio ITE hearing aids	Naída Lumity BTE hearing aids	Sky Lumity Pediatric BTE hearing aids
Form factor Tip: Battery size or rechargeability is indicated in the icons like this:	 Audéo I-Sphere Audéo I-R	 CROS I-R	 Virto I-Titanium Virto I-10 NW 0	 L-PR L-SP L-UP	 L-M L-PR L-SP L-UP
Charger	 Charger RIC Infinio ChargerGo RIC Infinio ChargerGo RIC Sphere Infinio	 Charger RIC Infinio ChargerGo RIC Infinio	Not applicable	 Phonak Charger Combi BTE 2 (Phonak Naída and Sky L-PR only)	 Phonak Pack (on-the-go charging)
Sound delivery system: Receivers Tip: See the receiver size and strength like this:	 SDS 6.0 receivers Sizes: 00, 0, 1, 2, 3 Strengths: S, M, P, UP (cShell)	 CROS Tube 6.0 Sizes: 0, 1, 2, 3	Not applicable		
Sound delivery system: Tubing and earhooks	Not applicable			 Slim Tube 4.0 or 5.0 Sizes: 00, 0, 1, 2, 3 or Standard tubing with these earhooks HE11 680 (damped) HE11 (undamped)	 Slim Tube 4.0 or 5.0 Sizes: 00, 0, 1, 2, 3 or Standard tubing with these earhooks HE11 680 (damped) HE11 (undamped) *Available as mini-hooks with additional colors.
Earpieces	 Domes Open, Vented, Power (Sizes: S, M, L) Cap (one size) SlimTip (acrylic, silicone) cShell	 Domes Open, Vented, Power (Sizes: S, M, L) CROS Tip SlimTip (acrylic, silicone)	Not applicable	 Domes Open, Vented, Power (Sizes: S, M, L) Cap (one size) SlimTip (acrylic, silicone) Ear mold (acrylic, silicone – different types possible)	
Wax filters	 CeruStop	Not applicable	 CeruStop HF3 or HF4 (Other options include extended receiver tube, wax spring, ext. receiver tube with wax spring)	Not applicable	
Accessories * Virto I Titanium, Virto I 10 NW 0, CROS I-R not compatible	 TV Connector			 PartnerMic	 Remote Control
Roger Remote mics * Virto I Titanium, Virto I 10 NW 0 only compatible with telecoil and neck loop; CROS I-R not compatible	 Roger On			 Roger Table Mic	 Roger Select™

Phonak hardware compatibility.

Loyal Phonak clients may still be wearing:

	Audéo Lumity RIC hearing aids	CROS Lumity Transmitter	Slim Lumity RIC hearing aid	Audéo Paradise RIC hearing aids	CROS Paradise Transmitter	Virto Paradise ITE hearing aids
Form factor Tip: Battery size or rechargeability is indicated in the icons like this: 	 L-312 L-R L-RT L-RL	 CROS L-R	 Slim L-R	 P-312 P-R P-RT P-13T	 P-R P-13	 P-Titanium P-312 NW O P-10- NW O P-312
Charger	 Phonak Charger Ease (L-R, L-RT)	 Phonak Charger (L-RL)	 Phonak Slim Charger	 Phonak Mini Charger (P-R, P-RT)	 Phonak Charger Case (P-R, P-RT)	 Phonak Power pack
Sound delivery system: Receivers Tip: See the receiver size and strength like this: 	 SDS 4.0 or 5.0 receivers Sizes: 00, 0, 1, 2, 3 Strengths: S, M, P, UP (cShell)	 CROS Tube 4.0 or 5.0 Sizes: 0, 1, 2, 3	 SDS 4.0 or 5.0 receivers Sizes: 00, 0, 1, 2, 3 Strengths: S, M, P, UP (cShell)	 CROS Tube 4.0 or 5.0 Sizes: 0, 1, 2, 3	 CROS Tube 4.0 or 5.0 Sizes: 0, 1, 2, 3	Not applicable
Earpieces	 Domes Open, Vented, Power (Sizes: S, M, L) Cap (one size) SlimTip (acrylic, silicone) cShell	 Domes Open, Vented, Power (Sizes: S, M, L) CROS Tip SlimTip (acrylic, silicone)	 Domes Open, Vented, Power (Sizes: S, M, L) Cap (one size) SlimTip (acrylic, silicone)	 Domes Open, Vented, Power (Sizes: S, M, L) Cap (one size) SlimTip (acrylic, silicone) cShell	 Domes Open, Vented, Power (Sizes: S, M, L) CROS Tip SlimTip (acrylic, silicone)	Not applicable
Wax filters	 CeruStop (SDS 5.0) CeruShield (SDS 4.0)	Not applicable	 CeruStop (SDS 5.0) CeruShield (SDS 4.0)	Not applicable	Not applicable	 CeruStop HF3 or HF4 (Other options include extended receiver tube, wax spring, ext. receiver tube with wax spring)
Accessories * Virto P-Titanium, Virto 10 NW O, Virto P-312 NW O and all CROS transmitters not compatible	 TV Connector		 PartnerMic	 Remote Control		
Roger Remote mics * Virto P-Titanium, Virto P-10 NW O, Virto P-312 NW O only compatible with telecoil and neck loop; CROS transmitters not compatible	 Roger On		 Roger Table Mic	 Roger Select™		