

# Phonak Audéo™ R Infinio.

Redefining extraordinary sound



Infinio Ultra



Phonak hearing aids are  
powered by advanced AI  
designed to help people hear  
better in the real world.

With Infinio Ultra, we are once again pushing the  
boundaries of what AI can achieve in hearing aids.  
**Discover how Phonak Audéo R Infinio can help  
you unlock a world of new possibilities.**



# Hear the world more clearly, every step of the way

Soft whispers, joyful laughter, and heartfelt conversations — Audéo R Infinio brings every sound to life. Our AI technology automatically adapts to your surroundings, providing exceptional sound quality in different listening situations.

- Fully automated and precise — boosting your speech understanding in noisy situations compared to the previous generation\*<sup>1</sup>
- From Jazz to HipHop, Classical to Electronic: whatever your music style, with our AI trained sound classification system you enjoy every beat



\* For adults with mild to severe hearing loss.

<sup>1</sup> Wright, A. et al. (2025). AutoSense OS 7.0 improves speech understanding with highly rated sound quality for challenging listening environments. Phonak Field Study News. Retrieved from [www.phonak.com/evidence](http://www.phonak.com/evidence).

# Effortless care, clear sound

Hearing aid care has never been easier, thanks to the new **EasyGuard™ domes**, which protect your hearing aids from earwax. EasyGuard features a unique sealed membrane that prevents wax from entering the hearing aid while allowing sound to pass through.

- Easy to clean<sup>2</sup>, so you can enjoy clear sound with less interruptions
- Protects the essential parts of your hearing aids from wax, helping them last longer<sup>2</sup>
- Eco-friendly and easy-to-use packaging



<sup>2</sup> Based on Phonak internal data. Contact [claims@sonova.com](mailto:claims@sonova.com) if you are interested in more information.

# Seamless connection. Minimal effort

Whether you're enjoying music, taking calls or watching a movie, we've got you covered with seamless connectivity to all your Bluetooth® devices.

- Enjoy one-step pairing to myPhonak and your Bluetooth® devices\*\*
- Pair with up to 8 Bluetooth® devices, and have 2 simultaneously connected
- Instantly shift between music and voice and enjoy uninterrupted hands-free calls and streaming across up to double the distance compared to our previous generations
- One click connection to Roger™ microphones, thanks to RogerDirect™



\*\* Smartphones and tablets

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Sonova AG is under license.

# Designed to withstand real life challenges

Experience class-leading precision with our expertly Swiss-engineered hearing solutions, designed to tackle life's challenges with confidence.

- Up to 31 hours of battery life\*\*\*, allowing you to engage in your favorite activities without interruption – all day long
- Your go-to solution for ultimate water and dust resistance



\*\*\* The maximum battery runtime is determined by battery condition, hearing loss, environmental sound classification, streaming usage and activated feature set.

# Meet Audéo R Infinio

## Colours

- Sand Beige
- Sandalwood
- Chestnut
- Champagne

- Silver Grey
- Graphite Grey
- Velvet Black



## 90 Premium

Speech in loud noise

Music

Speech in car

Comfort in noise

Speech in noise

Calm situation

Streaming speech

Streaming music

## 70 Advanced

Music

Speech in car

Comfort in noise

Speech in noise

Calm situation

Streaming speech

Streaming music



# Phonak chargers



**Phonak ChargerGo RIC Infinio**  
with built in battery for charging  
on the go



**Phonak Charger RIC Infinio**  
for cable charging

## 50 Standard

Speech in car

Comfort in noise

Speech in noise

Calm situation

Streaming speech

Streaming music

## 30 Essential

Speech in noise

Calm situation

Streaming speech

Streaming music

# Maximise your hearing potential with Roger

Functioning like a third ear, Roger picks up the voice of the speaker, television or telephone and wirelessly transmits it to your hearing aids, giving an extra boost in noisy environments and over distance.

With up to 65% improvement of word recognition in noise compared to using hearing aids alone<sup>3</sup>, Roger technology empowers you to communicate effectively in challenging situations.




<sup>3</sup> Zanin, J., Vaisberg, J., Swann, S., & Rance, G. (2025). Evaluating benefits of remote microphone technology for adults with hearing loss using behavioural and predictive metrics. *International Journal of Audiology*, 64(4), 1–9.

# myPhonak and myPhonak Junior

## Industry-leading apps

- Find My Hearing Aids, for peace of mind
- Custom programs, cleaning reminders, streaming balance, split volume control
- Ability to disable and enable streaming in the app without affecting Roger



# life is on

Since 1947, Phonak is dedicated to preserving social, emotional and physical life quality by opening new acoustic worlds. We believe that well-hearing equates to well-being and thus is essential for living life to the fullest. Today we offer the broadest portfolio of innovative hearing solutions. And, together with our hearing care professionals, we keep on focusing on what matters most: improving speech understanding, changing people's lives and having a positive effect on society as a whole.

[phonak.com](https://phonak.com)

028-2769-48/V1.00/2025-10/jb © 2025 Sonova AG. All rights reserved.

