

Reduce service appointments by up to 38% with EasyGuard¹.

Introducing Phonak EasyGuard™

Featuring a **sound-transmitting, sealed membrane** that protects the receiver from earwax.

Reduces the need for wax filter changes and helps to extend the lifespan of hearing aids, making care simple and convenient for your clients.

EasyGuard – Effortless care, clear sound



Comfortable silicone shell in multiple sizes
Made from the same material as our standard domes, for a comfortable fit

Wax-accumulating ribs
Can be easily wiped with a soft, dry cloth to remove accumulated wax

Sound transmitting membrane
Covers the receiver port to reduce the likelihood of wax or moisture damage

EasyGuard can help:



Reduce service appointments
by up to 38%



Free up time in your
clinic schedule



Increase client
satisfaction




Extend the lifespan
of hearing aids

In a 4 month home trial, 90% of participants indicated comfort levels matching or exceeding their traditional dome²


¹ Based on Phonak internal data. Contact claims@sonova.com if you are interested in more information.

² Sanchez, C. et Adler M.. (2025) New patented EasyGuard™ domes were preferred over traditional domes after a 4-month home trial. Phonak Field Study News. Retrieved from www.phonak.com/evidence

EasyGuard makes hearing aid care easier than ever



EasyGuard



Traditional dome
with wax filter

Comfortable for all
day wear



Sealed membrane for
additional protection
from earwax



1-step cleaning process
for ultimate convenience



Helps extend the lifespan
of your hearing aids



Changed every 3 months
or as needed



Dome only



Dome + wax filter