


Phonak hardware compatibility.

Our latest hearing aids

	Audéo Infinio RIC hearing aids		CROS Infinio Transmitter	Virto Infinio ITE hearing aids			Naída Lumity BTE hearing aids				Sky Lumity Pediatric BTE hearing aids			
Form factor Tip: Battery size or rechargeability is indicated in the icons like this: 	 Audéo I-Sphere Audéo I-R		 CROS I-R	 R Titanium 10 NW 0			 L-M L-PR L-SP L-UP				 L-M L-PR L-SP L-UP			
Performance level	70 and 90 30, 50, 70, 90		Not applicable	30, 50, 70, 90 70 and 90 30, 50, 70, 90			30, 50, 70, 90				50, 70, 90			
Charger	 Charger RIC Infinio ChargerGo RIC Infinio ChargerGo RIC Sphere Infinio		 Charger RIC Infinio ChargerGo RIC Infinio	 Charger Virto I-R			 Phonak Charger Combi BTE 2 (Phonak Naída and Sky L-PR only)				 Phonak Pack (on-the-go charging)			
Sound delivery system: Receivers Tip: See the receiver size and strength like this: 	 SDS 6.0 receivers Sizes: 00, 0, 1, 2, 3 Strengths: S, M, P, UP (cShell)		 CROS Tube 6.0 Sizes: 0, 1, 2, 3	Not applicable			Not applicable				Not applicable			
Sound delivery system: Tubing and earhooks	Not applicable			Not applicable			 Slim Tube 4.0 or 5.0 Sizes: 00, 0, 1, 2, 3 or Standard tubing with these earhooks HE11 680 (damped) HE11 (undamped)				 Slim Tube 4.0 or 5.0 Sizes: 00, 0, 1, 2, 3 or Standard tubing with these earhooks HE11 680 (damped) HE11 (undamped) <small>*Available as mini-hooks with additional colors.</small>			
Earpieces	 Domes Open, Vented, Power (Sizes: S, M, L) Cap (one size), EasyGuard™ (Open, Vented and Power) SlimTip (acrylic, silicone) cShell ¹		 Domes Open, Vented, Power (Sizes: S, M, L), EasyGuard™ (Open, Vented and Power) CROS Tip SlimTip (acrylic, silicone)	Not applicable			 Domes Open, Vented, Power (Sizes: S, M, L) Cap (one size) EasyGuard™ (Open, Vented and Power) SlimTip (acrylic, silicone) Ear mold (acrylic, silicone – different types possible)							
Wax filters	 CeruStop		Not applicable			 CeruStop HF3 or HF4 (Other options include extended receiver tube, wax spring, ext. receiver tube with wax spring)				Not applicable				
Accessories	* Virto I Titanium, Virto I 10 NW 0, CROS I-R not compatible		 TV Connector			 PartnerMic			 Remote Control					
Roger Remote mics	* Virto I Titanium, Virto I 10 NW 0 only compatible with telecoil and neck loop; CROS I-R not compatible		 Roger On			 Roger Table Mic			 Roger Select™					

¹ Length option 00 not available

Phonak hardware compatibility.

Loyal Phonak clients may still be wearing:

	Audéo Lumity RIC hearing aids	CROS Lumity Transmitter	Slim Lumity RIC hearing aid	Audéo Paradise RIC hearing aids	CROS Paradise Transmitter	Virto Paradise ITE hearing aids	
Form factor Tip: Battery size or rechargeability is indicated in the icons like this: 	 <p>L-312 L-R L-RT L-RL</p>	 <p>CROS L-R</p>	 <p>Slim L-R</p>	 <p>P-312 P-R P-RT P-13T</p>	 <p>P-R P-13</p>	 <p>P-Titanium P-312 NW 0 P-10- NW 0 P-312</p>	
Performance level	30, 50, 70, 90		30, 50, 70, 90	30, 50, 70, 90		30, 50, 70, 90	
Charger	 <p>Phonak Charger Ease (L-R, L-RT)</p>	 <p>Phonak Charger (L-RL)</p>	 <p>Phonak Slim Charger</p>	 <p>Phonak Mini Charger (P-R, P-RT)</p>	 <p>Phonak Charger Case (P-R, P-RT)</p>	 <p>Phonak Power pack</p>	
Sound delivery system: Receivers Tip: See the receiver size and strength like this: 	 <p>SDS 4.0 or 5.0 receivers Sizes: 0, 1, 2, 3 Strengths: S, M, P, UP (cShell)</p>	 <p>CROS Tube 4.0 or 5.0 Sizes: 0, 1, 2, 3</p>	 <p>SDS 4.0 or 5.0 receivers Sizes: 0, 1, 2, 3 Strengths: S, M, P, UP (cShell)</p>	 <p>CROS Tube 4.0 or 5.0 Sizes: 0, 1, 2, 3</p>	Not applicable		
Earpieces	 <p>Domes Open, Vented, Power (Sizes: S, M, L) Cap (one size)</p> <p>SlimTip (acrylic, silicone)</p> <p>cShell</p>	 <p>Domes Open, Vented, Power (Sizes: S, M, L)</p> <p>CROS Tip</p> <p>SlimTip (acrylic, silicone)</p>	 <p>Domes Open, Vented, Power (Sizes: S, M, L) Cap (one size)</p> <p>SlimTip (acrylic, silicone)</p>	 <p>Domes Open, Vented, Power (Sizes: S, M, L) Cap (one size)</p> <p>SlimTip (acrylic, silicone)</p> <p>cShell</p>	 <p>Domes Open, Vented, Power (Sizes: S, M, L)</p> <p>CROS Tip</p> <p>SlimTip (acrylic, silicone)</p>	Not applicable	
Wax filters	 <p>CeruStop (SDS 5.0) CeruShield (SDS 4.0)</p>	Not applicable	 <p>CeruStop (SDS 5.0) CeruShield (SDS 4.0)</p>	Not applicable		 <p>CeruStop HF3 or HF4 (Other options include extended receiver tube, wax spring, ext. receiver tube with wax spring)</p>	
Accessories * Virto P-Titanium, Virto 10 NW 0, Virto P-312 NW 0 and all CROS transmitters not compatible	 <p>TV Connector</p>		 <p>PartnerMic</p>	 <p>Remote Control</p>			
Roger Remote mics * Virto P-Titanium, Virto P-10 NW 0, Virto P-312 NW 0 only compatible with telecoil and neck loop; CROS transmitters not compatible	 <p>Roger On</p>		 <p>Roger Table Mic</p>	 <p>Roger Select™</p>			