

Phonak EON Feature Summary.

90 Premium		70 Advanced		50 Standard		30 Essential R	
Sphere	Automatic Spheric Speech Clarity 3.0	Automatic Spheric Speech Clarity 3.0	Automatic Spheric Speech Clarity 3.0	Spheric Speech Clarity 3.0	Manual Activation		
	Automatic Activation		Automatic Activation				
AutoSense OS™ AI 8.0	Highest Strength	High Strength	High Strength		Moderate Strength		
	Calm Situation		Calm Situation		Calm Situation		
	Spheric Speech in Loud Noise		Spheric Speech in Loud Noise		Speech in Noise		
	Speech in Car		Speech in Car		Speech in Car		
	Comfort in Noise		Comfort in Noise		Comfort in Noise		
	Music		Music		Streaming Speech		
	Streaming Speech		Streaming Speech		Streaming Music		
SmartSpeech Technology	Streaming Music	Streaming Music	Streaming Music	Motion Sensor Hearing	Spheric Speech Clarity 3.0	Roger Technology	Roger Technology
	SpeechSensor		SpeechSensor				
	Speech Enhancer		SpeechEnhancer				
	StereoZoom 2.0 (manual)		StereoZoom 2.0 (manual)				
	Dynamic Noise Cancellation		Dynamic Noise Cancellation				
	Motion Sensor Hearing		Motion Sensor Hearing				
	Spheric Speech Clarity 3.0		Spheric Speech Clarity 3.0				
	Roger Technology		Roger Technology				
	Speech in Loud Noise (manual)		Speech in Loud Noise (manual)				
	Speech in 360° (manual)		Speech in 360° (manual)				
	Tap Control		Tap Control				
	SoundRelax		SoundRelax				
	DuoPhone		DuoPhone				
	Real Ear Sound		Real Ear Sound				
	UltraZoom		UltraZoom				
	WindBlock 2.0		WindBlock 2.0				
	SoundRecover2		SoundRecover2				
	WhistleBlock		WhistleBlock				
	Tinnitus Balance		Tinnitus Balance				
	20 Channels		20 Channels		16 Channels		12 Channels



All levels: Universal Connectivity, now with Auracast™ broadcast streaming: Experience remarkable sound and handsfree calls from any Bluetooth® enable device, and experience the front row with Auracast™.

Phonak EON

Features and Benefits

WindBlock 2.0	WindBlock 2.0 reduces wind noise by up to 50 dB, improving the SNR by up to 30 dB resulting in improved target transparency, clearer speech and an overall more natural listening experience ¹
Adaptive Phonak Digital 4.0	A Phonak proprietary fitting formula ensuring a personalised first fit experience by enhancing overall sound quality based on individual hearing needs
AutoSense OS™ AI 8.0	Harnesses the power of machine learning to accurately identify and adjust to your sound environment – automatic when they need it, personalised through myPhonak where & when you want it
Spheric Speech Clarity 3.0	Instantly separates speech from noise, enabling you to enjoy clear speech from any direction, even in the most challenging group conversations ²
SpeechSensor	Improves hearing from the side and behind when in noise ³
Speech Enhancer	Makes it easier to understand soft speech or speech from a distance in quiet ⁴
Dynamic Noise Cancellation	Spatial noise cancellation to improve speech understanding in background noise ⁵
Motion Sensor Hearing	Automatically adjusts to optimise hearing while walking in noisy environments
Universal Connectivity	Universal Connectivity, now with Auracast™ broadcast streaming: Experience remarkable sound and handsfree calls from any Bluetooth® enable device, and experience the front row with Auracast™

What's the benefit of upgrading?

30 → 50	• Speech understanding in noisy environments from every direction with Spheric Speech Clarity (available as a manual programme) • Experience less background noise for a more comfortable experience when in car, plane or train • Automatically adjusts for better hearing when walking and talking
50 → 70	• Automatic activation of Spheric Speech Clarity, with increased denoising strength for even clearer speech in noisy environments from every direction • Dynamic noise cancellation to improve speech understanding in background noise • Dedicated automatic programme for live music
70 → 90	• Get the maximum denoising from Spheric Speech Clarity for the best speech understanding in noise from any direction from Phonak • Speech Enhancer makes it easier to understand soft speech in quiet, even from another room

¹ Vaisberg, J., Brunner, J., & Ulloa-Sanchez, D. (2026). Hearing confidently in wind: WindBlock 2.0 delivers comfort, natural sound quality and clearer speech. Phonak Field Study News retrieved from www.phonak.com/evidence

² Ulloa Sanchez et al. (2026). Phonak Audéo EON Sphere outperforms four other hearing aids for clear speech in loud noise. Phonak Insight retrieved from <https://www.phonak.com/evidence>

³ Latzel, M., Lesimple, C., & Woodward, J. (2022). New implementation of directional beamforming configurations show improved speech understanding and reduced listening effort. Phonak Field Study News retrieved from www.phonak.com/evidence.

⁴ Latzel, M., Lesimple, C., & Woodward, J. (2023). Speech Enhancer significantly reduces listening effort and increases intelligibility for speech from a distance. Phonak Field Study News retrieved from www.phonak.com/evidence.

⁵ Appleton, J. (2020). AutoSense OS 4.0 – significantly less listening effort and preferred for speech intelligibility. Phonak Field Study News retrieved from www.phonak.com/evidence.